


Cannulas

DETAILS OF Kelman Cystotome Cannula

Catalog No	AE-7101
Name	Kelman Cystotome Cannula
Category	Cannulas
Description	22 gauge, Overall length:32mm,Blunt on all sides, sharp at the tip

Image of Kelman Cystotome Cannula



DETAILS OF Lacrimal Cannula

Catalog No	AE-7010
Name	Lacrimal Cannula
Category	Cannulas
Description	23 gauge, straight 12mm malleable tip, Overall length:32mm

Image of Lacrimal Cannula

