


Cataract-Micro Coaxial Set

DETAILS OF Miyoshi Side port Manipulator

Catalog No	AE-2498
Name	Miyoshi Side port Manipulator
Category	Manipulators
Description	Features sandblasting on the ball and shaft for better control in the presence of viscoelastic during the Counter Traction Implantation Technique(See AE-9036SP) • Ball tip diameter: 0.35mm • Angle to tip: 1.65mm

Image of Miyoshi Side port Manipulator



DETAILS OF Akahoshi Capsulorrhexis Cross-Action Forceps

Catalog No	AE-4345
Name	Akahoshi Capsulorrhexis Cross-Action Forceps
Category	Forceps
Description	Overall length: 113mm, Tip length: 12mm, Tip angled at 135 degrees. Unique cross action mechanism reduces loss of viscoelastic material through the incision. Consecutively deep anterior chamber can be maintained and rhexis edge can be better controlled. The pointed tip can make the initial puncture of the capsule easily and grasps the capsule firmly. Can also be used through a side port incision.

Image of Akahoshi Capsulorrhexis Cross-Action Forceps

