


Conductive Keratoplasty

DETAILS OF Milne CK Light Touch Marker points @ 6, 7, and 8

Catalog No	AE-2845
Name	Milne CK Light Touch Marker points @ 6, 7, and 8
Category	Markers
Description	Overall length: 103mm, Outer Diameter: 8.5mm, Inner Diameter: 5mm. Unique extension design allows re-marking of cornea in offset manner of 22.5 degrees to original CK treatment marks (when performing a CK enhancement). Features a built-in sturdy centering device for precise placement on the cornea, offering excellent visibility of the cornea compared to other traditional barrel markers. Features a built-in sturdy centering device for precise placement on the cornea, offering excellent visibility of the cornea compared to other traditional barrel markers. 24 point, with 8 points each at 6, 7 and 8mm

Image of Milne CK Light Touch Marker points @ 6, 7, and 8



DETAILS OF Milne CK marker 16 points @ 7 and 8

Catalog No	AE-2847
Name	Milne CK marker 16 points @ 7 and 8
Category	Markers
Description	Overall length: 103mm, Outer Diameter: 8.5mm, Inner Diameter: 5mm. Unique extension design allows re-marking of cornea in offset manner of 22.5 degrees to original CK treatment marks (when performing a CK enhancement). Features a built-in sturdy centering device for precise placement on the cornea, offering excellent visibility of the cornea compared to other traditional barrel markers. Features a built-in sturdy centering device for precise placement on the cornea, offering excellent visibility of the cornea compared to other traditional barrel markers. 16 point, with 8 points each at 7 and 8mm

Image of Milne CK marker 16 points @ 7 and 8

