

Corneal Transplant

DETAILS OF John DXEK / DSAEK Fixation Hook

Catalog No	AE-2182
Name	John DXEK / DSAEK Fixation Hook
Category	Hooks
Description	Exquisite design allows for effective fixation of the donor "taco" during air injection to unfold the donor disk within the recipient anterior chamber (AC)

Image of John DXEK / DSAEK Fixation Hook



DETAILS OF John DXEK / DSAEK Descmet Stripper

Catalog No	AE-2874
Name	John DXEK / DSAEK Descmet Stripper
Category	Spatulas
Description	Overall length: 110mm, Tip length: 16mm. Unique design allows for excellent apposition to the inner-dome curvature of the recipient cornea, allowing for uniform detachment of the host DM without damage to the corneal stroma. Curvilinear design is effective in reaching the DM over all inner curvatures of the recipient cornea. Round handle permits rotational movement along the long axis of the instrument.

Image of John DXEK / DSAEK Descmet Stripper

