

Glaucoma

DETAILS OF Castroviejo Caliper - Straight

Catalog No	AE-1500
Name	Castroviejo Caliper - Straight
Category	Calipers
Description	Straight; Overall length:86mm, Measures 1.0mm-20.0mm in 1.0mm increments

Image of Castroviejo Caliper - Straight



DETAILS OF Crandall USC Glaucoma Scleral Planer

Catalog No	AE-2430
Name	Crandall USC Glaucoma Scleral Planer
Category	Glaucoma Planer
Description	Overall length:103mm, Angle to tip:10mm, Tip width:1.00mm, Dual cuttings edges. Rapidly remove layers of sclera to unroof schlemm's canal for viscocanalostomy and other non-invasive trabeculotomy procedures

Image of Crandall USC Glaucoma Scleral Planer

