


# Intralase

## DETAILS OF Lehmann Intralase Spatula

Catalog No	AE-2827
Name	Lehmann Intralase Spatula
Category	Spatulas
Description	Overall length: 113mm, Shaft length: 10mm, Tip bend: 135 degrees, Tip length: 11mm. The smooth tip enables easy entry under the flap without adjacent epithelial damage or disruption. Minimal technique learning curve and greatly eases surgeon's ability to lift the flap while keeping bed dry. While delicate compared to other instruments for this purpose, it is superior to cannulas both in regards to sturdiness and sterility. Spatula works well to find and mark or begin entry into the flap margin for enhancements.

Image of Lehmann Intralase Spatula



## DETAILS OF Lehmann Intralase Spatula

Catalog No	AE-2827S
Name	Lehmann Intralase Spatula
Category	Spatulas
Description	Lehmann Intralase Spatula with semi sharp edge. Overall length: 113mm, Shaft length: 10mm, Tip bend: 135 degrees, Tip length: 11mm. The smooth tip enables easy entry under the flap without adjacent epithelial damage or disruption. Minimal technique learning curve and greatly eases surgeon's ability to lift the flap while keeping bed dry. While delicate compared to other instruments for this purpose, it is superior to cannulas both in regards to sturdiness and sterility. Spatula works well to find and mark or begin entry into the flap margin for enhancements.

Image of Lehmann Intralase Spatula

