

LRI (Limbal Relaxing Incision)

DETAILS OF Arnold Guarded Speculum

Catalog No	AE-1060
Name	Arnold Guarded Speculum
Category	Specula
Description	The blades of the speculum are "guarded" to seclude 90% of the lid margin and are smoothly finished to allow relatively easy insertion for greater patient comfort. The blades are non-reflective and the screw mechanism prevents lid squeezing when opened to the full potential of the patient's palpebral fissure.

Image of Arnold Guarded Speculum



DETAILS OF Zaldivar LRI Axis Marker

Catalog No	AE-2862
Name	Zaldivar LRI Axis Marker
Category	Markers
Description	Overall length: 105mm. Created to mark the cornea at 10.5mm. To avoid error due to the torsional movement of the globe, fixed marks are located on the 90-degree meridian as reference points for marking in the upright position and referencing in the recumbent position. Corrects astigmatism of one and two diopters. Designed with an axis extension to easily align with the steepest meridian while working with a degree gauge for precise alignment.

Image of Zaldivar LRI Axis Marker

