


# Lasik

## DETAILS OF Bansal LASIK Marker

Catalog No	AE-2800
Name	Bansal LASIK Marker
Category	Markers
Description	Overall length:105mm, Noncontact OZ: 6.0mm, Radials:3, Pararadials:2

Image of Bansal LASIK Marker



## DETAILS OF McPherson Tying Forceps

Catalog No	AE-4373
Name	McPherson Tying Forceps
Category	Forceps
Description	Overall length:102mm. Straight long handle, Smooth jaws, 5mm long platform. Ideal for 8-0 to 11-0 sutures.

Image of McPherson Tying Forceps

