


Phakic IOLs

DETAILS OF IOL Manipulator and Iris Retractor

Catalog No	AE-2218
Name	IOL Manipulator and Iris Retractor
Category	Manipulators
Description	Overall length:115mm, Tip width:.80mm. Ideal for side to side and up and down lens manipulation and to stretch anterior capsule

Image of IOL Manipulator and Iris Retractor



DETAILS OF Clearpath Diamond Knife

Catalog No	AE-8158
Name	Clearpath Diamond Knife
Category	Diamond Knives
Description	Angled Handle 2.5mm blade,

Image of Clearpath Diamond Knife

