


Retina

DETAILS OF VR-TIS 23 Gauge Spatula

Catalog No	AE-2917
Name	VR-TIS 23 Gauge Spatula
Category	Retina
Description	Overall length: 152mm, Tip length: 40mm. Flat, thin curved spatula, 2.75mm tip length designed for blunt dissection and delamination of membrane from the retina. Gentle curve allows manipulation along the surface of the retina.

Image of VR-TIS 23 Gauge Spatula



DETAILS OF VR-TIS 23 Gauge Curved Scissors

Catalog No	AE-5606
Name	VR-TIS 23 Gauge Curved Scissors
Category	Retina
Description	23 gauge. Overall length: 146mm, Tip length: 35mm. Multi purpose scissors featuring 1.25mm curved blades. This allows the blade to follow the profile of the globe for cutting the membrane on the surface .

Image of VR-TIS 23 Gauge Curved Scissors

