


Toric IOL

DETAILS OF Nuijts Toric Axis Mkr w/5mm ctr opening

Catalog No	AE-2740
Name	Nuijts Toric Axis Mkr w/5mm ctr opening
Category	Markers
Description	Used with a Mendez Degree Gauge to accurately mark the steepest meridian • Circular opening maximizes visibility during marking • Laser lines on the non-marking surface indicate marking points • Lightweight, all titanium body

Image of Nuijts Toric Axis Mkr w/5mm ctr opening



DETAILS OF Pre-op Toric Reference Marker

Catalog No	AE-2793D
Name	Pre-op Toric Reference Marker
Category	Markers
Description	Marking pattern is located on both sides for marking the left and right eye temporally. Lightweight, all titanium body.

Image of Pre-op Toric Reference Marker

