

# ColorMouse-RM™

# Bierley

**Non-Glare  
Image**



**Full Colour**

**White on  
Black**

**Black on  
White**

The ColorMouse-RM is ideally for reading or looking at magazines, maps, prescription bottles, food labels, stamps, photos, and much more. By simply plug the ColorMouse-RM into any regular television, connect the power and press the large blue button, your television will change to the correct channel automatically on European models. For the USA and rest of the world just press your remote control "Video" button. The RM (Reverse Mode) function makes it easier on the eyes for those people with acute light sensitivity.

- Compatible with any TV using video input
- High contrast image, less glare
- Vivid full colour mode
- Ergonomic and easy to use
- High contrast negative image (pure white text on black background)
- High contrast positive image (pure black text on white background)
- 24x Magnification on a 20" TV :(High)
- 13x Magnification on a 20" TV :(Low)
- SCART connection with Auto-Detect (Europe)
- RCA connection (USA & rest of world)
- 6 ~ 9 feet of easy store cable
- Durable, CE approved, built to last

**1  
Year  
Warranty**

[www.bierley.com](http://www.bierley.com)

**CALL NOW**

**0800 0430282**