

ColorMouse™

Bierley

**Vivid
Colour**



Full Colour

- See anything in true full colour
- Ergonomic and easy to use

The ColorMouse is ideal for reading books or looking at magazines, maps, prescription bottles, food labels, stamps, photos and much more. Simply plug the ColorMouse into any regular television, switch it on and enjoy reading even the smallest text in large, sharp, vivid colour. Operation couldn't be easier because your television will automatically switch to the correct channel for using your ColorMouse and will automatically switch back to the channel that you were watching once you switch the ColorMouse off. This feature is unique to our European SCART products.

You can leave your ColorMouse connected to the television so that it is always ready to use.

- Compatible with any TV using video input
- 24x Magnification on a 20" TV :(High)
- 13x Magnification on a 20" TV :(Low)
- SCART connection with Auto-Detect (Europe)
- RCA connection (USA & rest of world)
- 6 ~ 9 feet of easy store cable
- Durable, CE approved, built to last

**1
Year
Warranty**

www.bierley.com

CALL NOW

0800 0430282