

MonoMouse-RM™

Bierley

**Non-Glare
Image**



**White on
Black**

**Black on
White**

The MonoMouse-RM provides three modes for reading, namely regular mode, reverse mode (pure white text on black background), and enhanced mode (pure black text on white background). With quick installation by plugging the MonoMouse-RM into any regular television and pressing the large blue button, the magnifier is ready to be used and your television becomes a powerful low vision aid. Operation couldn't be easier because your television will automatically switch to the correct channel for using your MonoMouse-RM and will automatically switch back to the channel that you were watching once you switch the MonoMouse-RM off. This feature is unique to our European SCART products.

You can leave your MonoMouse-RM connected to the television so that it is always ready to use.

- Compatible with any TV using video input
- High contrast image, less glare
- Ergonomic and easy to use
- 24x Magnification on a 20" TV :(High)
- High contrast negative image (pure white text on black background)
- High contrast positive image (pure black text on white background)
- 24x Magnification on a 20" TV :(High)
- 13x Magnification on a 20" TV :(Low)
- SCART connection with Auto-Detect (Europe)
- RCA connection (USA & rest of world)
- 6 ~ 9 feet of easy store cable
- Durable, CE approved, built to last

**1
Year
Warranty**

www.bierley.com

CALL NOW

0800 0430282