




# Optical Products 2008



PO Box 611 Christchurch NEW ZEALAND  
Tel: +64 3 366 2999 Fax: 365 2072  
E-mail: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)  
[www.blaxall.net](http://www.blaxall.net)

**blaxall optics ltd**

# Close Focus Monoculars



*Leading Edge Optics*

Very compact and lightweight. All models have anti reflection coated optics with a spiral twist focus. They will focus from close to infinity. A rubber grip gives ease of use. Complete with soft protective case and straps.



M1025D

**BLAXALL OPTICS LTD**

P0 Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

ITEM	POWER	LENS SIZE	CLOSE FOCUS	FOV
M320	3x	20mm	20cm	11.5°
M412	4x	12mm	19cm	12.5°
M616	6x	16mm	26cm	9.5°
M821	8x	21mm	40cm	7.2°
M1025D	10x	25mm	35cm	5.2°

# Clip-on Monoculars



*Leading Edge Optics*

**BLAXALL OPTICS LTD**

P0 Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

Monoculars designed to clip onto spectacle frames. Diopter adjustable twist focus. Model M2525C is achromatic with an anodized aluminum body. The advantage of these monoculars is that they are hands free and relatively inexpensive. They are also beneficial for prolonged distance viewing and have a wide field of view.



ITEM	POWER	LENS SIZE	TYPE	FOV
M2525C	2.5x	25mm	Achromatic	11.5°
M2828C	2.8x	28mm	Standard	12.5°

# Deluxe ALLBLAX LED MAGNIFIERS



MG6979



BLAXALL OPTICS LTD

PO Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

The latest generation of innovative LED magnifiers providing exceptional long battery life and a bright image with good contrast.

The aspheric lenses are coated to improve light transmission. Designed for the partially sighted with ease of use an objective. The heads may be purchased separately and are interchangeable.

## TORCH

ITEM	MAGNIFICATION	DIOPTER	SIZE
MG6951	4x	12D	70mm
MG6952	7x	23D	50mm
MG6953	8x	28D	35mm
MG6956	10x	36D	35mm
MG6958	5x	20D	60mm
MG6955	6x	21D	57mm
MG6957	3x	8D	100x75mm
MG6959	3x	8D	100x50mm

## POCKET

ITEM	MAGNIFICATION	DIOPTER	SIZE
MG6970	4x	12D	70mm
MG6971	7x	23D	50mm
MG6976	6x	21D	57mm
MG6972	8x	28D	35mm
MG6973	10x	36D	35mm
MG6974	12x	44D	35mm
MG6975	14x	50D	35mm
MG6978	5x	20D	60mm
MG6977	3x	8D	50x75mm
MG6979	3.5x	10D	50x75mm

# Hybrid LED TORCH MAGNIFIERS



This range differs from the standard range in the aspect that the LED diode is twice as bright giving even greater contrast with no loss of diode life. The hybrid LED handle can be purchased separately and is fully interchangeable on all heads.



**BLAXALL OPTICS LTD**

PO Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

ITEM	MAGNIFICATION	DIOPTR	SIZE
MG6951D	4x	12D	70mm
MG6952D	7x	23D	50mm
MG6953D	8x	28D	35mm
MG6956D	10x	36D	35mm
MG6958D	5x	20D	60mm
MG6957D	3x	8D	100x75mm
MG6959D	3x	8D	100x50mm

# ALLBLAX “OYSTER” MAGNIFIERS



*Leading Edge Optics*

**BLAXALL OPTICS LTD**

PO Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

Folding pocket magnifiers for inspection, field work,  
low vision and daily use. Complete with pouch.



ITEM	POWER	DIOPTR	LENS SIZE	TYPE
<b>MG7171</b>	4x	12D	35mm	Aspheric
<b>MG7172</b>	6x	20D	35mm	Aspheric
<b>MG7173</b>	8x	28D	35mm	Aspheric
<b>MG7174</b>	10x	39D	30mm	Aspheric


# ALLBLAX™

## HIGH QUALITY MAGNIFIERS



*Leading Edge Optics*

A new hybrid range of lightweight acrylic magnifiers with exceptional optics. Crafted with advanced technology, incorporating sophisticated plastics. Designed for ease of handling and as the frames are completely transparent, the orientation field is not impaired. Developed with the assistance of The Lighthouse, New York.



**BLAXALL OPTICS LTD**

PO Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

ITEM	POWER	LENS SIZE	WEIGHT
MGV50	4x	50 mm	40 gm
MGV60	3x	60 mm	45 gm
MGV80	2.5x	80 mm	75 gm
MGV100	2x	100 mm	110 gm
MGV100R	2x	100/50 mm	100 gm

# Allblax™ High Power Magnifiers



*Leading Edge Optics*

**BLAXALL OPTICS LTD**

PO Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

Lightweight aspheric hand held magnifiers.  
Affordable magnification ideal for the majority  
of low vision patients plus a range of other  
applications.



ITEM	POWER	LENS SIZE	DIOPTR
MGV100HP	3x	100x75mm	8D
MGV70HP	4x	70mm	12D
MGV60HP	5x	60mm	16D
MGV57HP	6x	57mm	20D
MGV50HP	7x	50mm	24D

# LED RIMLESS MAGNIFIER



*Leading Edge Optics*

High quality acrylic magnifiers with ergonomic rubber covered handle. Rimless magnifiers provide distraction free viewing. The two LED diodes provide exceptionally bright illumination. Complete with batteries and soft pouch. Blister packed with counter display stand available.



**BLAXALL OPTICS LTD**

P0 Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

ITEM	POWER	LENS SIZE	ILLUMINATION
<b>MGRL85</b>	2x,6x	85x80mm	LED
<b>MGRL75</b>	3x	75mm	LED
<b>MGRL90</b>	2x	90mm	LED
<b>MGRL11</b>	2x	110mm	LED

# BRIGHT FIELD MAGNIFIERS




*Leading Edge Optics*

Bright Field (dome) magnifiers provide distortion-free viewing and have a calculated light guidance that directs all possible illumination onto the object. This provides a brighter image without an additional light source. Useful to use in combination with other optical aids.



MG6912



MG6913

**BLAXALL OPTICS LTD**

PO Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

ITEM	POWER	DIOPTR	SIZE
MG6913	3x	8D	60mm
MG6914	3x	8D	75mm
MG6912	4x	12D	65mm
MG6911	4x	12D	80mm
MG6908	3x	10D	95mm

# LOW VISION SPORTS GLASSES



*Leading Edge Optics*

Galilean low vision sports glasses for the visually impaired. Excellent coated glass optics. Each eye can be focused individually and the nose bridge is adjustable. Each unit is supplied with case and side eye shields.



**BLAXALL OPTICS LTD**

PO Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

ITEM	POWER	DIOPTR	LENS SIZE
B92828	2.8x	7 D	28mm
B92525	2.5x	5.1 D	25mm (achro)
B93528	3.5x	10 D	25mm (achro)

# BAR MAGNIFIERS



*Leading Edge Optics*

Acrylic Bar magnifiers with a guide line for easy guiding. These magnifiers allow strain-free reading over long periods. Complete with protective case.



MG1008

**BLAXALL OPTICS LTD**

PO Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

ITEM	MAGNIFICATION	SIZE
MG1008	2.5x	125x25mm
MG1008L	2.5x/5x	210x25mm

# MISCELLANEOUS



*Leading Edge Optics*

Stand magnifiers especially useful for older people who may have hand trembles. The stand on item MG75100 has a folding stand.



MG75100



**BLAXALL OPTICS LTD**

PO Box 611  
Christchurch

Phone: +64 3 366 2999

Fax: +64 3 365 2072

Email: [blaxall@paradise.net.nz](mailto:blaxall@paradise.net.nz)

ITEM	MAGNIFICATION	SIZE
MG75100	3x	75x100mm
MG835H	8x	35mm
MG750H	7x	50mm

