

BinoSpecs

SeeTV

Hands-Free Glasses for Intermediate and Distance Viewing

Maximize the size of TV screens with the easiest-to-fit, binocular, telescopic distance system ever introduced - SeeTV glasses! These hands-free, head-mounted glasses provide 2.1x magnification.

SeeTV glasses will focus on objects from distances of 10 feet to infinity, and are perfect for watching TV, sporting events, movies, theatre, bird watching, or any other distance viewing activity in which a magnified image would be helpful.

With flexible temples, the SeeTV telescopic glasses are the most comfortable television viewing device ever. Don't spend thousands on a large screen TV - try the SeeTV first!

Features of SeeTV glasses:

Galilean design offers a bright, wide field of view (20 degrees) and 2.1x magnification;

PXM optics and flexible plastic temples allow for light weight (only 1.7 ounces) with an extremely comfortable fit;

Attractive, modern, hi-tech frame;

Frosted objective lens housing eliminates distracting reflections;

Moulded bridge, with no sharp edges, makes the system comfortable to wear, even for long periods of time;

Engineered for wearers with binocular PDs of 60mm to 68mm;

Tech notes:

Comes with soft case which doubles as a cleaning cloth

Magnification: 2.1x

Lens Size: 32mm

Lens Type: Galilean

Lens Material: PXM

Temple Length: 135mm

Bridge: Saddle

Frame Size: 137mm

Focal Range: 10 feet to infinity

Visual Field: 20 degrees

PD Range: 60mm - 68mm

Weight: 1.7 oz



blaxall optics ltd

Tel: +64 3 366 2999 Fax: 03 365 2072 Email: blaxall@paradise.net.nz