

Accessories

LENS DRILLING UNIT



With notch cutting function

Drill work frequency : 6000 – 13000 rev/min

Power supply : 100-120V/60Hz AC or 200-240V/50Hz AC

Acceptable Lens Materials : CR-39, polycarbonate, glass (mineral)

Acceptable processing diameter : 0.8 mm – 2.8 mm

ELECTRICAL SOLDERING UNIT



Electrical soldering unit

Power supply : 110 or 220 v

Average consumption : 220 w

Dimensions : 20 x 25 x 15 cm

Weight 10 Kg.

ACCESSORIES FOR SOLDERING UNIT



Spare Carbon for
soldering unit



Flux with 25g
liquid



Soldering string 10
m length



Pistol for professional
soldering unit

ELECTRICAL SET FOR FRAME FINISHING OR RE-TOUCHING



Power supply Bivolt

Dimensions : 15,5 cm

Weight: 80 gr.

LENS COLOURING MACHINE



Lens colouring machine with 6 pots

Machine for lens colouring and treatment.

Constant temperature controlled by a micro-processor.

High performance for uniform lens colouring and treatment.

Power supply 110 or 220 v

Capacity: 8 litres

Average consumption 1300 w

Dimensions : 27 x 39 x 27 cm

Weight: 8,5 Kg.

BURTON'S LAMP



Burton's lamp

Composite system with special bifocal lens (enlarges vision).

Power supply 110 or 220 v

Dimensions : 16,5 x 23 x 6 cm

Weight 850 gr.

Lamps 2 white + 2 ultraviolet

Visor square or round

DOUBLE SIDED MIRROR



Table mirror with or without enlarging side.

Dimensions : 37 x 27 x 22 cm

Diameter .22,5 cm

Weight 950 gr.

Colours : red, white, light blue,

silver, beige, brown, black

OVAL FIRENZE MIRROR



Oval Firenze mirror

Oval mirror with elegant finishing without enlarging side.


Dimensions : 42 x 22 x 17 cm

Diameter 22 x 32 cm

Weight 2 Kg.

Colours : steel - onyx

TRIAL LENS SET



Complete set of trial lenses with metal rim and aluminium case.

Lens diameter : 36,5 mm


PLIERS




001: Square Pad Adjusting Pliers




002: Round Pad Adjusting Pliers.



003: Pliers for removal of non "clip on" type nose pads.



004: Holding Pliers to settle glasses bows in acetate frames.




005: Holding Pliers to settle glasses bows in metal frames.



006 : Thin Pointed Pliers to adjust the internal side of the frame.



007 : Thin Pointed Pliers for general and pad arm adjustment.



008 : Flat Nose Pliers for bridges and joints adjustment.



009 : Round Nose Pliers for general and pad arm adjustment.



010 : Curved Nose Pliers to shape and bend delicate components or not easily reachable parts.



011 : Side Bended Pliers for frames bows adjustment and inclination



012 : Size Testing Pliers for testing lens size during finishing process



013 : Lens Aligning Pliers to adjust the differences in lens axis after mounting them in the frame



014 : Pliers for frame bending for a better lens setting.



015 : Side Cutting Pliers to cut screws, lenses and general adjustments. Small or big size



016 : Front Cutting Pliers specific for screws



017 : Pliers for cristal lenses chipping Small or big size

TOOL SET



Tool set

Complete set of pliers, screwdrivers and tools in an elegant brief-case.

SCREWDRIVER AND BLADES SET



Screwdriver and blades set

Complete with 10 different blades
(Phillips, screwdriver and nut driver)
Rigid case for storage

TWEEZERS



THIN POINTED : N. 1, 2, 3, 4 or 5

POINTED : to grip screws and nuts

SOLDER TWEEZER

HOLD TWEEZER

EXTRACTOR TWEEZER

JEWEL 2V



Buffing Unit Jewel 2V
Lens and frame polishing unit

Features :

useful for cleaning and polishing lens and frame. Complete with a rag buffing wheel, felt wheel and rough stick. Preset for the aspirator.

2 speeds :

1400- 2800 rpm

Power : 275 W

Dimensions :

Shaft length : 440 mm

Height : 200 mm

Depth : 140 mm

Weight :

9 kg

In accordance with EC norms.

COMET



Features :

small dimension frame heater with electronic temperature regulation. Aluminium structure. Complete with all accessories (cup for glasses, diffusion device).

Power supply :

220 V, 50 Hz- 110 V, 60 Hz

Dimensions :

29x13x10 cm

Weight :

2 kg

In accordance with EC norms.

STAR VENTILETTA



Features :

Frame heater with electronic temperature regulation.

Power supply :

100-120V/60Hz AC or 200-240V/50Hz AC

Power : 1500W

Temperature : 140~150° C

Dimensions :

31x25x16 cm

Weight :

2 kg

In accordance with EC norms.

