


AO-Type Projector Slides

LEA Symbols®, LEA Numbers®, Sloan Letters



Designed for use in all manual Chart Projectors including products from R.H. Burton, Marco, Reichert, Topcon, and Woodylyn.

A 3200VAC presents LEA Symbols® in 19 lines from 20/125 to 20/10 with multiple lines from 20/50 to 20/10 and Sloan Letters in 20 lines from 20/400 to 20/10 with multiple lines from 20/50 to 20/10.

3200VAC

B 3201VAC presents LEA Numbers® in 17 lines from 20/125 to 20/10 with multiple lines from 20/50 to 20/10 and Sloan Letters in 20 lines from 20/400 to 20/10 with multiple lines from 20/50 to 20/10.

3201VAC

C 3202VAC presents LEA Symbols® (18 lines), LEA Numbers® and Sloan Letters (16 lines) from 20/100 to 20/10 with multiple lines from 20/50 to 20/10.

3202VAC

D 3203VAC presents LEA Symbols® in 20 lines from 20/200 to 20/10 with multiple lines from 20/50 to 20/10 and Sloan Letters in 20 lines from 20/250 to 20/10 with multiple lines from 20/50 to 20/10.

3203VAC

E 3204VAC presents LEA Symbols® in 11 lines from 20/200 to 20/16 and Sloan Letters in 17 lines from 20/400 to 20/16 with multiple lines from 20/40 to 20/20.

3204VAC

F 3205VAC presents LEA Symbols®, LEA Numbers® (10 lines), and Sloan Letters (14 lines) from 20/125 to 20/16. Multiple lines from 20/40 to 20/20 for Sloan Letters.

3205VAC

G 3206VAC presents LEA Numbers® in 11 lines from 20/200 to 20/16 and Sloan Letters in 17 lines from 20/400 to 20/16 with multiple lines from 20/40 to 20/20.

3206VAC

H 1077VAC Pediatric presents Allen figures from 20/200 to 20/30, Tumbling E from 20/200 to 20/20, and Landolt C from 20/200 to 20/20.

1077VAC

I 1103VAC Astigmatism & Vergence presents Tumbling E from 20/400 to 20/200, Sloan Letters from 20/100 to 20/10, Astigmatic Clock Test, Vertical Lines for Vergence from 20/40 to 20/10, Sloan Letters from 20/50 to 20/10, and Astigmatic Blocks.

1103VAC

J 1113VAC Acuity Astigmatism presents Tumbling E and Sloan Letters from 20/300 to 20/10, Astigmatic Clock Test, Sloan Numbers from 20/50 to 20/10, and Tumbling E from 20/200 to 20/15.

1113VAC

K 1180VAC Adult Acuity presents combination Sloan Letters and Sloan Numbers from 20/400 to 20/10, Sloan Numbers from 20/70 to 20/10, and Tumbling E from 20/200 to 20/15.

1180VAC

L 1192VAC Adult/Pediatric Acuity presents combination Sloan Letters and Sloan Numbers from 20/400 to 20/15, Tumbling E from 20/200 to 20/25, and Allen Figures from 20/200 to 20/30.

1192VAC

M 1193VAC General Acuity presents Tumbling E and Sloan Letters from 20/400 to 20/10, Astigmatic Clock Test, Sloan Letters from 20/50 to 20/10, Tumbling E from 20/50 to 20/15, and Allen Figures from 20/100 to 20/40.

1193VAC

N 1234VAC Adult Vectorgraphic presents Sloan Letters Monocular from 20/100 to 20/15, Astigmatic Clock, Acuity Balance from 20/30 to 20/20, Acuity Suppression from 20/40 to 20/25, Binocular Acuity from 20/30 to 20/15, Fixation Disparity (with/without fusion lock), and Stereopsis with 260-73 seconds of arc.

1234VAC