


# Color Vision Tests

Color Vision Books



## Pseudoisochromatic Plate Ishihara Compatible (PIPIC) Color Vision Test 24-Plate Edition

by Terrace Waggoner, OD

Quickly screens for normal color vision and both congenital and acquired color vision defects. Classifies to type (protan, deutan, and tritan). The only test available that contains both tritan (blue/yellow deficiency) and basic color vision test plates. Provides additional plates to comply with testing requirements of the military and most other government agencies.

730009 .....



## Color Vision Testing Made Easy

by Terrace Waggoner, OD

This test uses simple, easy-to-recognize symbols and is appropriate for children as young as age 3 years. It is also the gold standard for testing individuals with developmental challenges. As a bonus feature, Quick 6 Color Test is included for a limited time—absolutely FREE!

730001 .....

Available in three different versions!



## Richmond HRR Standard Pseudoisochromatic Test

Reengineered by Jay Neitz, PhD, and James Bailey, OD, PhD

Quickly screens for normal color vision and congenital and acquired color vision defects, classifying to type—protan/deutan (red/green) and tritan (blue/yellow)—and extent of defect. Reengineering and clinical results reviewed by the International Color Vision Society. Results submitted for publication in Visual Neuroscience. Structure identical to original AO HRR.

730006 .....



## Quick 6 Color Test

Abbreviated quick red/green color deficiency screening test.

730008 .....

## Richmond Laminated HRR Test

This test is ideal for pediatric testing, because it uses simple shapes (O, X, Δ) that are easily recognized by children. Special nonglare laminate keeps colors from becoming distorted and prevents fingerprint smudges. Reversible occluding glasses are included.

730005 .....

## Richmond HRR/Amsler Combination Laminated Color Book

Testing for acquired color vision defects is now easier than ever with two of the key tests combined in one book. The new HRR 4th edition is combined with a book of Amsler Charts, including Diagonals, Red Lines, and Subdivided Squares. The set also includes reversible occluding glasses.

13396LAR .....



## Improved Good-Lite® Pseudoisochromatic Color Vision Test

This improved test offers 8 numbers and 8 pediatric symbols with larger testing plates for easy screening of children and adults—even those with low vision. Checks for both red/green and blue/yellow deficiencies. Includes a basic color vision test plate, user-friendly instructions, tabbed-sections, and protective and antiglare coating on each test page.

730000 .....