



Fixation Products

Attention-grabbing Items




Small and Medium Fixation Paddles

The small 2" (5 cm) diameter "smiling face" fixation paddle was developed based on studies by Robert Fanz. Features include grating on reverse side for preferential looking, small pictures of a bunny and a giraffe on end of stick for fixation, and marks for measuring pupil size and corneal diameter measurement on the handle.

The medium 5" (13 cm) diameter fixation paddle is for assessing children with severe disabilities to determine the distance at which the children respond to a high-contrast face figure.

- A 253100 Medium
- B 253000 Small



Toucan Fixation PD Rule

Toucan fixation and PD rule with near optotypes, including LEA Symbols®, LEA Numbers®, and HOTV Letters for 20/32, 20/60, 20/80, and 20/100. Includes inch/mm rules and diopter scale plus pupillometer in half mm increments. 6" long.

4970R



Optokinetic Flag

This flag is used to produce optokinetic nystagmus (OKN) physiologic jerk in a patient's eyes. This is elicited when the

patient attempts to fixate on moving, alternative, contrasting squares that are moved in the field of vision. Comes with instructions. 37.25" (94.62 cm) long.

4856R



Distance Fixation Animals

With animation, noise, and lights, these distance Fixation Animals hold any child's attention. They also come attached to a shelf for easy mounting. A new remote control foot pedal is an option to the traditional wired foot pedal. You can also choose from a handheld remote control or an animal without a pedal or remote. Both remote controls have multiple channels.

Frog Drummer (10-1/2"H x 9"W) is loud, raucous, and attention-grabbing. The time-honored Husky and Sheep are best suited where the fixation distance is less than 15 feet. The Husky and Retriever Puppy are most easily recognized by younger children. The Retriever Puppy can be used as a replacement animal for COW industries units. Advise us if you will have multiple animals in one office.

Animal Type	Remote Control	Remote Pedal	Wired Pedal	No Pedal, No Remote
A Rolly Poly Cow	181CR	181RPR	181R	181NFR
B Retriever Puppy	136CR	136RPR	136R	136NFR
C Duck	174CR	174RPR	174R	174NFR
D Piglet	175CR	175RPR	175R	175NFR
E Elephant	173CR	173RPR	173R	173NFR
F Sheep	190CR	190RPR	190R	190NFR
G Frog Drummer	193CR	193RPR	193R	193NFR
H Siberian Husky	178CR	178RPR	178R	178NFR
Price Eash				



Flashing Fixation Star

For mounting on a retinoscope or ophthalmoscope; includes 3 flashing star fixations

with 9 mm magnetic disks in box and double-sided tape for mounting. Fits most retinoscopes.

4890R



Near Point Fixation Animals

These fun, near point fixations are both colorful and easily recognized by

children. Pin lever in the back permits the introduction of motion. Heat-resistant nylon washer allows figures to fit over a 3/16" transilluminator (muscle light) neck. These are intended for use at 16" to 20". Set of 2 animals.

4960R