

# HOTV Visual Acuity

Distance and Near Visual Acuity Products

## Distance Wall Charts



### HOTV Folding Chart (10 Feet, 3 Meters)

For screening/testing the vision of children and adults ages 3 years and older. Proportionally-spaced lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent. 50% spaced rectangle (Massachusetts Visual Acuity Test format) with line sizes from 20/100 to 20/8 (6/30 to 6/2.4) equivalent. Has a hole for hanging on the wall. Includes response key, flash cards, and instructions. Chart measures 18" x 18" (45.7 cm x 45.7 cm) and folds to 18" x 9" (45.7 cm x 22.9 cm).


600302 .....



### HOTV Chart (10 Feet, 3 Meters)

For screening/testing the vision of children and adults 3 years and older. Proportionally-spaced lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent. 50% spaced rectangle (Massachusetts Visual Acuity Test format) with line sizes from 20/100 to 20/8 (6/30 to 6/2.4) equivalent. Has a hole for hanging on the wall. Includes response key, flash cards, and instructions. Chart measures 18" x 18" (45.7 cm x 45.7 cm).

600303 .....



### Combination LEA Symbols®/ HOTV MassTest™ Wall Chart (10 Feet, 3 Meters)

This test is designed specifically for the Massachusetts Preschool Vision Screening Protocol. LEA Symbols® are on one side; HOTV are on the opposite side. One practice line, two line sizes of 20/40 (equivalent) for 3 year olds, and two line sizes of 20/32 (equivalent) for 4+ year olds. Has a hole for hanging on the wall. Includes chart, response key, flash cards, and a copy of the Massachusetts Vision Screening Protocol. Chart measures 10" x 10" (25.4 cm x 25.4 cm).

259700 .....

## Handheld Flip Charts



### MassTest™ Protocol Flip Chart (10 Feet, 3 Meters)

This test is designed specifically for the Massachusetts Preschool Vision Screening Protocol. One practice card, two line sizes of 20/40 for 3 year olds, and two line sizes of 20/32 for 4+ year olds. Includes response key, flash cards, and a copy of the Massachusetts Vision Screening Protocol. Pages are offset for easy flipping. 4 pages, 3.25" x 9" (8.25 cm x 22.9 cm).

600304 .....



### Preschool SiGHT® LEA Symbols® and HOTV Screening System

Specially developed to screen for anisometric amblyopia in children ages 3 through 5. Uses both LEA Symbols® and HOTV letters in line sizes 20/100, 20/50, and 20/30 equivalent. 40" (1 meter) cord with attached occluder ensures proper testing distance. Includes a set of children's fun horse occluders, complete instructions, recording form, response key, and flash cards. 13 pages. 4.25" x 4.5" (10.8 cm x 11.43 cm).

A 400001 Crowded .....

B 400002 Uncrowded .....



### Single HOTV Letter Book With 50% Spaced Crowding Bars (10 Feet, 3 Meters)

This book features 50% contour interaction bars to provide single optotype acuity with the proper crowding effect. Optotype lines sizes from 20/100 to 20/16 (6/30 to 6/4.8) equivalent. Includes response panel, flash cards, and instructions. 17 pages. 4" x 4" (10.2 x 10.2 cm).

600350 .....