

HOTV Visual Acuity

Distance and Near Visual Acuity Products

Professional Cabinet Charts

Chart Size: 10' x 18", 25.4 cm x 45.7 cm



Crowded HOTV Translucent Chart (10 Feet, 3 Meters)

Crowded lines; line sizes range from 20/160 to 20/16 equivalent. Set includes chart, response key, flash cards, and instructions.

800004 Set

800726 Chart Only



HOTV Translucent Chart (20 Feet, 6 Meters)

Uncrowded HOTV chart with line sizes ranging from 20/100 to 20/16. Response key and flash cards are not included.

800704

Near Vision Tests



MassTest™ Near Vision Protocol Screener With 16-inch (40 cm) Cord

This test is designed specifically for the Massachusetts screening protocol. 50% spaced rectangle (Massachusetts Visual Acuity Test format) with line sizes 20/100 and 20/32 equivalent. HOTV printed on one side and Sloan Letters on the opposite side. The response key is printed on the test card. 6" x 10".

740000



HOTV Near Vision Card With 16-inch (40 cm) Cord

This double-sided card combines proportionally-spaced lines on one side and 50% spaced rectangle (Massachusetts Visual Acuity Test Format) on the other. Includes response key, flash cards, and instructions. 8" x 10" (20.3 cm x 25.4 cm).

731000



Crowded HOTV Near Vision Card

Crowded line sizes from 20/125 to 20/20 (6/38 to 6/6) equivalent. Instructions

printed on reverse. Includes two paper isolating masks. 5.5" x 8" (14 cm x 20.3 cm).

726000

Replacement Components




A HOTV Response Key

700523

B HOTV Flash Cards

700515



Combination LEA Symbols® and HOTV Response Key

700030



Combination LEA Symbols® and HOTV Flash Cards

700031