



Illuminated Cabinets

Distance Visual Acuity Products



ESV1200 Quantum™ Cabinet*

LED lights provide uniform non-calibrated illumination across the entire testing surface. Automatically adjusts to incoming power of 100 to 240 VAC.¹ Glare lights and wireless remote control are not provided with this unit. Can be powered by an optional battery pack. Charts are not provided with this unit and must be ordered separately. 9.25"W x 15"H x 3.25"D (38.1 cm x 23.5 cm x 8.3 cm). CE and CSA Listed.

ESV1200

ESV1500 Quantum™ Cabinet*

Using advanced patented technology, this cabinet self-calibrates illumination on the chart surface for both photopic (85 cdm²) and mesopic (3 cdm²) light level testing. LED lights provide uniform and highly calibrated illumination across the entire testing surface (Patent 5,078,486). Automatically adjusts to incoming power of 100 to 240 VAC.¹ Controlled by infrared wireless remote. Optional floor stand and glare lights (see page 31) are also available. 4.5 lbs. Can be powered by an optional battery pack. Charts are not provided with this unit and must be ordered separately. 9.25"W x 15"H x 3.25"D (38.1 cm x 23.5 cm x 8.3 cm). CE and CSA Listed.

ESV1500

* Compatible with all existing 9" x 14" (23 cm x 33.5 cm) charts and stands.

¹ Power plug adaptor may be needed for outlets outside the United States.



Model A+ Translucent Eye Cabinet*

The Model A+ Translucent Eye Cabinet combines all of the best features of the previous Models A and M, portability, and storage. With an added hinged-lid storage compartment (for extra charts), this versatile cabinet can be placed on a table, hung on a wall, or set on a stand to accommodate most screening or testing situations. A wall mount bracket is also available. Housed in a lightweight (5.8 lbs.) metal cabinet. The model A+ will provide years of accurate and dependable service. Charts are not provided with this unit and must be ordered separately. 9"W x 14"H x 4.25"D (22.9 cm x 35.6 cm x 10.8 cm). CSA Listed.

600600 (120VAC, 60Hz, 0.18A)

Replacement Bulb

White Fluorescent Replacement Bulb

Used with Model A+, Model A, or Model M. F8T5D 8 watt

600500



Model A Translucent Cabinet*

The lightweight (5 lbs.), versatile Model A can be placed on a table, hung on a wall, or set on a stand to accommodate accurate 10' to 20' visual acuity screening and testing. Housed in a durable metal cabinet, the Model A provides dependable service and has been the most popular cabinet of practitioners for more than 25 years. A wall mount bracket is also available. Uses F8T5D 8-watt fluorescent bulb. Charts are not provided with this unit and must be ordered separately. 9"W x 15"H x 4"D (22.9 cm x 38.1 cm x 10.2 cm). CSA Listed.

600000 (120VAC, 60Hz, 0.18A)

Isolating Masks



Isolating Masks

Set of two with 1" x 1" and 1.25" x 1.25" holes. Fits ESV1500, ESV1200, Model A+, or A.

600076



Linear Test Mask

For ESV1500, ESV1200, Model A+, or A.

600014