

The LEA Test System includes line tests, crowded symbol tests, and single symbol tests to allow assessment of the visual system across these three functional areas. Over the years, the number of LEA tests has increased to meet the screening and assessment needs of children and adults at different ages and developmental levels.

Distance Wall Charts



MassTest™ Folding Distance Chart (10 Feet, 3 Meters)

For screening or testing individuals beginning at age 3. 50% spaced rectangle (Massachusetts Visual Acuity Test format) with line sizes from 20/200 to 20/8 (6/60 to 6/2.4) equivalent. Has hole for hanging the chart on the wall. Includes response key, flash cards, and instructions. 18" x 20" (45.7 cm x 50.8 cm). Folds to 18" x 10" (45.7 cm x 25.4 cm).

258800




Make vision screening fun!
See page 34 for pediatric occluders.



15-Line Distance Chart (10 Feet, 3 Meters)

For screening/testing children ages 4 years and older. Proportionally-spaced lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent. 50% spaced rectangle (Massachusetts Visual Acuity Test format) with line sizes from 20/100 to 20/8 (6/30 to 6/2.4) equivalent. Has a hole for hanging the chart on the wall. Includes response key, flash cards, and instructions. 17" x 21" (43 cm x 53.4 cm).

257000



10-Line Distance Chart (10 Feet, 3 Meters)

For screening children beginning at age 3. Proportionally-spaced lines; line sizes range from 20/200 to 20/16 (6/60 to 6/4.8) equivalent. The 20/160 (6/48) and 20/100 (6/30) lines have been omitted. 50% spaced rectangle (Massachusetts Visual Acuity Test format) with line sizes from 20/80 to 20/16 (6/24 to 6/4.8) equivalent. Has a hole for hanging the chart on the wall. Includes response key, flash cards, and instructions. 16" x 16" (40.6 cm x 40.6 cm).

256000

10' and 20' notations only (not shown)

256011



15-Line Folding Distance Chart (10 Feet, 3 Meters)

For screening or testing individuals beginning at age 3. Proportionally-spaced lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent. 50% spaced rectangle (Massachusetts Visual Acuity Test format) with line sizes from 20/100 to 20/8 (6/30 to 6/2.4) equivalent. Has a hole for hanging the chart on the wall. Includes response key, flash cards, and instructions. 17" x 21" (43 cm x 53.4 cm). Folds to 17" x 10.5" (43 cm x 26.7 cm).

250100



10-Line Folding Distance Chart (10 Feet, 3 Meters)

For screening children beginning at age 3. Proportionally-spaced lines; line sizes range from 20/200 to 20/16 (6/60 to 6/4.8) equivalent. The 20/160 (6/48) and 20/100 (6/30) lines have been omitted. 50% spaced rectangle (Massachusetts Visual Acuity Test format) with line sizes from 20/80 to 20/16 (6/24 to 6/4.8) equivalent. Has a hole for hanging the chart on the wall. Includes response key, flash cards, and instructions. 16" x 16" (40.6 cm x 40.6 cm). Folds to 16" x 8" (40.6 cm x 20.3 cm).

250200

10' and 20' notations only (not shown)

250211