



### Low Vision LEA Symbols® Chart Set (40 inches, 1 Meter)

This low vision test is designed for examining children and adults with severe visual impairments. The largest symbols are 50M in size, corresponding to 1/50 (20/1000, 6/300, 0.02). The smallest symbols are 1M in size, corresponding 1/1 (20/20, 6/6, 1.0) at a distance of one meter. The test comes with attached cord to maintain the proper testing distance. Includes response key, flash cards, and instructions. Designed courtesy of August Colenbrander, MD. Chart measures 16" x 21.5" (40.7 cm x 54.6 cm). Folds to 16" x 10.75" (40.7 cm x 27.3 cm).


258000 .....



### Low Vision LEA Numbers® Chart (40 inches, 1 Meter)

This low vision test is designed for examining children and adults with severe visual impairments. The largest numbers are 50M in size corresponding to 1/50 (20/1000, 6/300, 0.02) and the smallest symbols are 1M in size corresponding 1/1 (20/20, 6/6, 1.0) at a distance of one meter. The test comes with an attached cord to maintain the proper testing distance. Includes instructions. Designed courtesy of August Colenbrander, MD. 16" x 21.5" (40.7 cm x 54.6 cm). Folds to 16" x 10.75" (40.7 cm x 27.3 cm).

272100 .....



### Runge Low Vision Near Card (16 inch, 40 cm)

Sloan Letters with line sizes ranging from 20/1000 to 20/16 (equivalent). Near acuity measurements and their equivalents are listed on reverse side. 13" x 7.5".

756900 .....



### Lea Symbols® Low Vision Book

This low vision test is designed for examining children and adults with severe visual impairments. This complete book includes testing for 10' (3 meters), 5' (1.5 meters) and for near vision at 16" (40 cm). Distance testing line sizes range from 20/1600 to 20/16 equivalent. Near vision testing line sizes range from 20/800 to 20/16 equivalent. This book comes bound either on the left side or top. Butterfly occluder glasses for children and metal frame occluder glasses for adults are included.

5274R Left-Side Bound .....

5275R Top Bound .....



### Lea Numbers® Low Vision Book

This low vision test is designed for examining children and adults with severe visual impairments. This complete book includes testing for 10' (3 meters), 5' (1.5 meters) and 16" (40 cm) near vision. Distance testing line sizes range from 20/1600 to 20/16 equivalent. Near vision testing line sizes range from 20/800 to 20/16 equivalent. This book comes bound either on the left side or top. Butterfly occluder glasses for children and metal frame occluder glasses for adults are included.

5131R Left-Side Bound .....

5273R Top Bound .....