



Reading Visual Acuity

Distance and Near Vision Acuity Products

Continuous Text and Reading Cards



Lebensohn Near Vision Card (14 inch, 35.6 cm)

This 2-sided contemporary reading card, designed by James Lebensohn, MD, provides a full range of near vision exercises, including words, numbers, and OX symbols. Optotypes include 20/80, 20/400, 20/260, 20/200, 20/100, 20/65, 20/50, 20/40, 20/30, 20/25, and 20/20. Notation equivalents include decimal readings in English and metric, Jaeger, point size, and visual efficiency percentages. Red/green bars employed on 5 pt to 12 pt optotypes. Card size is 5.5" x 8.5" (14 cm x 22 cm).

11901R.....



Contemporary Near Vision Card (16 inch, 40 cm)

This 2-sided contemporary reading card provides familiar reading exercises, such as a crossword puzzle, help wanted, football and basketball scores, music song bar, calendar, playing cards, recipes, instructions, telephone white page listings, and TV listings. Continuous format reading examples provided for 20/160, 20/100, 20/80, 20/50, 20/40 and 20/30 equivalent. Notation equivalents include Jaeger and point size. Card size is 8.5" x 11" (21.5 cm x 28 cm).

11900R.....



Rosenbaum Numbers Near Vision Pocket Card (14 inch, 35.6 cm)

The Rosenbaum Near Vision Numbers Pocket Card includes all the traditional Rosenbaum features plus near point lines in HOTV, LEA Numbers®, and LEA Symbols® for 20/100, 20/90, 20/80, and 20/32. Also features an Amsler grid (black and white lines), an accommodation test, and pupilometer figures in both full and half circles to the edge for measuring pupils from 2 mm to 9 mm. Card fits in the pocket and shows point, Jaeger, and distance equivalents from 20/800 to 20/20. Uses numbers, directional "E", and "XO" patterns. 3.5" x 6.5" (9.5 cm x 15.9 cm).

5032R.....



Runge Pocket Near Vision Card (16 inch, 40 cm)

Sloan Letters with line sizes ranging from 20/500 to 20/16 (equivalent). Near acuity measurements and their equivalents are listed on the reverse side. Also available as LEA Symbols® (pg. 13) or LEA Numbers® (pg. 15). 3.75" x 6.25" (9.5 cm x 15.9 cm).

756000 Sloan Letters.....



Sloan Letters Near Vision Pocket Card (16 inch, 40 cm)

Sloan letters with line sizes of 20/400 to 20/10 equivalent on one side. Continuous text is printed on the reverse side with line sizes 20/80 to 20/16. Card measures 3.75" x 6.25" (9.5 cm x 16 cm).

744000.....



Fonda-Anderson Near Vision Reading Card (16 inch, 40 cm)

One side displays continuous text while the other side has LEA Symbols®, LEA Numbers®, and Sloan optotypes. Card measures 5" x 8" (12.7 cm x 20.3 cm).

727000.....



Dr. Good's Near Vision Reading Card (14 inch, 35.6 cm)

First-grade words. Line sizes from 20/160 to 20/20 (6/48 to 6/6) equivalent. Card measures 5.5" x 8" (14 cm x 20.3 cm).

724000.....