

# Sloan Letters Visual Acuity

Distance and Near Visual Acuity Products



## Distance Wall Chart



## Sloan Letters Folding Chart (10 Feet, 3 Meters)

Proportionally-spaced lines; line sizes range from 20/200 to 20/8 (6/60 to 6/2.4) equivalent. 18" x 18" (45.7 cm x 45.7 cm); folds to 18" x 9" (45.7 cm x 22.9 cm).

735000 .....



## Magnetic Sloan Letters Chart (10 Feet, 3 Meters)

This durable chart is made of plastic coated galvanized steel with 2 sets of 20 magnetic strips (40 total) with 5 characters on each strip. Meets all qualifying standards for military aviators (e.g., officer candidates, pilots, air crewmen, etc). Approved by the Bureau of Medicine (BUMED). Chart measures 9.5" x 11" (24.1 cm x 28 cm).

400500 20/20 size .....

400900 20/40 size (shown) .....

## ETDRS Translucent Distance and Low Contrast Charts

Chart Size: 24.5" x 25.5", 62 cm x 65 cm



## Original Series ETDRS Translucent Distance Charts (13 Feet, 4 Meters)

Proportionally-spaced (geometric progression) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent. Standardized format used in National Eye Institute clinical studies and in low vision evaluations that require accurate and repeatable visual acuity measurements.

500013 Chart "R" .....

500014 Chart "1" .....

500015 Chart "2" .....



## ETDRS Translucent Low Contrast Charts (13 Feet, 4 Meters)

Proportionally-spaced (geometric progression) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent.

500023 25% .....

500019 10% .....

500020 5% .....

500021 2.5% .....

500022 1.25% .....



## "2000" Series ETDRS Translucent Distance Charts (13 Feet, 4 Meters)

Proportionally-spaced (geometric progression) lines; line sizes range from 20/200 to 20/10 (6/60 to 6/3) equivalent. Revised version of the original ETDRS series.

500016 Chart "1" .....

500017 Chart "2" .....

500018 Chart "3" .....



## "2000" Series ETDRS Translucent Distance Charts (10 Feet, 3 Meters)

Proportionally-spaced (geometric progression) lines; line range sizes from 20/200 to 20/10 (6/60 to 6/3) equivalent. Revised version of the original ETDRS series.


500010 Chart "1" .....

500011 Chart "2" .....

500012 Chart "3" .....

## Professional Cabinet Charts

Chart Size: 10" x 18", 25.4 cm x 45.7 cm



## C Sloan Letters Chart (10 Feet, 3 Meters)

Linear-spaced lines; line sizes range from 20/160 to 20/16 (6/48 to 6/4.8) equivalent.

800713 .....

## D Sloan Letters Translucent Chart (10 Feet, 3 Meters)

Proportionally-spaced lines; line sizes range from 20/160 to 20/8 (6/48 to 6/2.4) equivalent.

800735 .....

**PLEASE NOTE:** All charts may also be mounted or hung on the wall.