

Tumbling E Visual Acuity

Distance and Near Visual Acuity Products



Quantum Series, Model A+, and Model A Translucent Charts

Chart Size: 9' x 14", 22.9 cm x 35.8 cm



Tumbling E Chart (10 Feet, 3 Meters)

Linear-spaced lines; line sizes range from 20/100 to 20/16 equivalent. Lines are labeled with 10' and 20' notations only.

600738



Tumbling E Chart (20 Feet, 6 Meters)

Wide-spaced lines; line sizes range from 20/100 to 20/16. Lines are labeled with 20' notations only.

600728



Tumbling E Chart (10 Feet, 3 Meters)

Wide-spaced lines; line sizes range from 20/100 to 20/16 equivalent. Lines are labeled with 10' and 20' notations only.

This chart is recommended by the Texas Department of Public Health for vision testing and screening.

600722



Tumbling E Chart (20 Feet, 6 Meters)

Linear-spaced lines; line sizes range from 20/80 to 20/16. Lines are labeled with 20' notations only.

600735

Professional Cabinet Charts

Chart Size: 10' x 18", 25.4 cm x 45.7 cm



Tumbling E Chart (10 Feet, 3 Meters)

Linear-spaced lines; line sizes range from 20/160 to 20/16 (6/48 to 6/4.8) equivalent.

800714



Tumbling E Chart (20 Feet, 6 Meters)

Linear-spaced lines; line sizes range from 20/100 to 20/20 (6/30 to 6/6).

800729