



H u m a n I n f o r m a t i o n M a n a g e m e n t S e r v i c e



#### Technical Specification

LCD	4.3" wide type
Magnification	3x ~ 13x(close-up camera)
Image storing	Capture and storing up to 20pcs
Color mode	Full color, black/white, white/black, yellow/blue, yellow/black, blue/white
Brightness	4 levels
Battery	Lithium Polymer (rechargeable) 4.5 hours continuous use
Size	17.8cm(W) × 7.8cm(D) × 2.8cm(H) 7"(W) × 3"(D) × 1.1"(H)
Weight	261g, 0.57lb



Fingertip Communication

URL <http://www.himsintl.com>



## SenseView Duo

SenseView Duo combines distance camera and close-up & handwriting camera in one device.



H u m a n I n f o r m a t i o n M a n a g e m e n t S e r v i c e



Fingertip Communication

# Compact, stylish and slim-sized portable CCTV



**Clear contrast image**  
Clear and high contrast image



**Easy handwriting**  
You can easily use handwriting function without additional devices.



### Distance Viewing

View objects in the distance such as street signs or classroom chalk boards.



### Dual camera

- Close-up camera : for viewing and handwriting mode
- Distance camera : for objects 1m (3.3feet) or further away



### Image storing

- Storing capacity : 20 images



### Various magnification rate

- Min. magnification : 3x
- Max. magnification : 13x



### Convenient screen interface

- Menu and panning/selection buttons for data management and battery level icon



### Image Panning

- Captured and stored images can be magnified and panned for clearer view.