



HOWARD INSTRUMENTS INC.



IrisMate[®] Iris Retractors

Flexible Translimbal Iris Retractors



Each Sterilization Case contains 6 Semi-disposable Iris Retractors

Affordable

Easy Insertion and Removal

All MEDICAL GRADE materials

Flexible nylon hook retracts iris with minimal damage to the tissue

Specially designed silicone cinch provides superior grip

Works through 0.3mm - 0.5mm self sealing incision

Safely used in Aphakic, Pseudophakic and Phakic Eyes