



# INFUSION CANNULAS

## DISPOSABLE INFUSION CANNULA

- 20 Gauge
- 45 degree beveled tip
- Luer lock adapter
- Uses 260mm of Silastic® tubing
- Length 270mm, 10.6 inches

### IS-1025

- 2.5mm

### IS-1040

- 4.0mm

### IS-1060

- 6.0mm



## METAL INFUSION CANNULAS

### HIGH VISCOSITY DISPOSABLE METAL INFUSION CANNULA

- 20 Gauge
- 45 degree beveled tip
- Can be used as standard infusion cannula
- Designed for easy injection of high viscosity substance such as silicone oil
- Tubing length 381mm, 15 inches

### MIS-1040

- 4.0mm tip

### MIS-1060

- 6.0mm tip



## SILICONE OIL / PERFLOURDECALINE

### SILICONE OIL

NOTFDAAPPROVED

#### SO-1000

- Sterile
- 1000 centistoke
- 10ml

#### SO-2000

- Sterile
- 2000 centistoke
- 10ml

#### SO-5000

- Sterile
- 5000 centistoke
- 10ml



### PERFLOURDECALINE

NOTFDAAPPROVED

#### PFD-1000

- Sterile
- 5ml

#### PFD-2000

- Sterile
- 7ml


Howard Instruments, Inc.  
 4749 Appletree  
 Tuscaloosa, AL 35405  
 Website: [howardinstruments.com](http://howardinstruments.com)  
 email: [howard@howardinstruments.com](mailto:howard@howardinstruments.com)  
 Phone: (205) 553-4453 Fax: (205) 556 9267



# Howard Instruments, Inc. *SPECIAL VITRO-RETINAL ACCESSORIES*



H-2800 Viscous Fluid Injector



H-2860 Perflourdecaline Injection  
Needle

## *SILICONE OIL / PERFLOURDECALINE*

### **SILICONE OIL** NOT FDA APPROVED

#### **SO-1000**

- Sterile
- 1000 centistoke
- 10ml

#### **SO-2000**

- 2000 centistoke

#### **SO-5000**

- 5000 centistoke



### **PERFLOURDECALINE** NOT FDA APPROVED

#### **PFD-1000**

- Sterile
- 5ml

#### **PFD-2000**

- 7ml




## *REUSABLE SCLERAL PLUGS & TABLET*



**R20-10000**  
20 gauge



**R19-10000**  
19 gauge





**R18-10000**  
18 gauge

Supplied with 4 plugs  
Color-coded  
Steam Autoclavable silicone rubber tablet  
Raised support for easy grip with forceps (not included)

## *DISPOSABLE SCLERAL PLUGS & TABLET*

### **RDP-1000**

Supplied with 4 plugs,  
2, 19 gauge  
2, 20 gauge  
Polypropylene tablet



### **RCA-1000** Scleral Plug Forceps

- Lightweight
- Developed for use with virtually any traditional plug set available
- Manufactured from the highest quality American stainless steel