



# Howard Instruments, Inc.

## *POLLY'S VITREOUS CUTTERS*



### **J-225-O20GA**

- 20 gauge
- 20psi
- For Ocutome I and II
- Port is 17/1000" from the tip
- Port size is adjustable
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided

### **J-225-O20GBN**

- 20 gauge
- 20 psi
- 90 degree bent needle
- For Ocutome type machines

### **J-225-O23G**

- 23 gauge
- 20 psi
- For Ocutome type machines
- Port is 17/1000" from the tip
- Double lumen tubing provided

### **J-225-O23GBN**

- 23 gauge
- 20 psi
- 90 degree bent needle
- For Ocutome type machines

### **J-225-O23GH**

- 23 gauge
- 20 psi
- For Ocutome type machines. Port is 17/1000" from the tip.
- Port is 17/1000" from the tip
- Hook at end of port for easier guiding of scar tissue into the cutting port
- Placement of hook often makes forceps unnessecary
- Used during the more delicate stages of scar tissue removal
- Used for primary debulking of scar tissue
- Double lumen tubing provided

### **J-225-O25G**

- 25 gauge
- 20 psi
- For Ocutome type machines
- Port is 17/1000" from the tip
- Double lumen tubing provided

### **J-225-10000**

- 20 gauge
- 20 psi
- For Alcon 10000/20000 type machines
- Port is 17/1000" from the tip
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided

### **J-225-AMO**

- 20 gauge
- 20 psi
- For AMO type machines
- Port is 17/1000" from the tip.
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided

### **J-225-C20G**

- 20 gauge
- 22 psi
- For Chiron type machines
- Port is 17/1000" from the tip
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided

### **J-225-INTRA**

- 20 gauge
- 30 psi
- For Intraoptic type machines
- Port is 17/1000" from the tip
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided

### **J-225-OMS**

- 20 gauge
- 12 psi
- For OMS type machines
- Port is 17/1000" from the tip
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided

Howard Instruments, Inc.  
4749 Appletree  
Tuscaloosa, AL 35405  
Website: [earshowardinstruments.com](http://earshowardinstruments.com)  
email: [howard@howardinstruments.com](mailto:howard@howardinstruments.com)  
Phone: (205) 553-4453 Fax: (205) 556 9267



# Howard Instruments, Inc.

## *POLLY'S VITREOUS CUTTERS*



### **J-225-SD**

- 20 gauge
- For Surgical Design type machines
- Port is 17/1000" from the tip
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided

### **J-225-I**

#### **Infusion Sleeve**

- For Anterior Vitrector

### **J-225-M20G**

- 20 gauge
- 30 psi
- For Microvit type machines, also for Daisy, Premier and Nidek type machines
- Port is 17/1000" from the tip
- Port size 0.9mm
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided



### **J-225-M20GA**

- 20 gauge
- 30psi
- For Microvit type machines, also for Daisy, Premier and Nidek type machines
- Port is 17/1000" from the tip
- Port size is adjustable
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided

### **J-225-M20GBN**

- 20 gauge
- 30 psi
- 90 degree bent needle
- For Microvit type machines

### **J-225-M23G**

- 23 gauge
- 30 psi
- For Microvit, Daisy and Premiere type machines
- Port is 17/1000" from the tip
- Double lumen tubing provided

### **J-225-M23GBN**

- 23 gauge
- 30 psi
- 90 degree bent needle
- For Microvit type machines

### **J-225-M23GH**

- 23 gauge
- 30 psi
- For Microvit, Daisy and Premiere type machines
- Port is 17/1000" from the tip
- Hook at end of port for easier guiding of scar tissue into the cutting port
- Placement of hook often makes forceps unnecessary
- Used during the more delicate stages of scar tissue removal
- Used for primary debulking of scar tissue
- Double lumen tubing provided

### **J-225-M25G**

- 25 gauge
- 30 psi
- For Microvit, Daisy and Premiere type machines
- Port is 17/1000" from the tip
- Double lumen tubing provided

### **J-225-O20G**

- 20 gauge
- 20 psi
- For Ocutome I and II
- Port 17/1000" from the tip
- Port size 0.9mm
- J-225-I, irrigation sleeve, needed for anterior segment surgery
- Double lumen tubing provided

Howard Instruments, Inc.  
4749 Appletree  
Tuscaloosa, AL 35405  
Website: [earshowardinstruments.com](http://earshowardinstruments.com)  
email: [howard@howardinstruments.com](mailto:howard@howardinstruments.com)  
Phone: (205) 553-4453 Fax: (205) 556 9267



Howard Instruments, Inc.

## *DISPOSABLE VITREOUS CUTTERS*



**J-225-DM**

- 30 psi
- For the Microvit and Premiere
- May also be used for DORC machine

**J-225-DMG**

- 30 psi
- For the Fritz Ruck machine

**J-225-DO**

- 20 psi
- For Ocutome 800 and 8000 machines



**J-225-DAM**

- 30 psi
- For Microvit and Premiere
- With irrigation sleeve

**J-225-PJD**

- 30 psi
- For PhacoJack machine
- With irrigation sleeve

**J-225-DA10000**

- 20 psi
- For Alcon 10000 and Alcon 20000 machine.
- With irrigation sleeve

## *POLLY'S VITREOUS MICROSCISSORS*

### *Vertical Microscissors*

**J-225-VS30**

- 30 psi
- For Microvit and similar machines

**J-225-VS20**

- For Ocutome and similar machines
- 20 psi

### *Horizontal Microscissors*

**J-225-HS30**

- For Microvit and similar machines
- 30 psi


**J-225-HS20**

- For Ocutome and similar machines
- 20 psi



**J-225-M20GBN**

- 20 gauge
- 30 psi
- 90 degree bent needle
- For Microvit type machines



Howard Instruments, Inc.  
4749 Appletree  
Tuscaloosa, AL 35405  
Website: [earshowardinstruments.com](http://earshowardinstruments.com)  
email: [howard@howardinstruments.com](mailto:howard@howardinstruments.com)  
Phone: (205) 553-4453 Fax: (205) 556 9267