




RETINAL ACCESSORIES



H-2115 Membrane Spatula
* Hooked
* 2.5 mm
* Separating fibrovascular membranes from the retina



H-2120 Membrane Spatula
* Curved
* All purpose membrane manipulation



H-2100 Membrane Peeler
* Blunt tip
* Locating and stripping membranes from the retinal surface



H-2710 Scleral Depressor
* Double ended
* 2.4 diameter conical tip
* Slightly curved 32 mm shaft
* 6.5 mm wide bar tin heavily curved 35 mm shaft



H-2135 Membrane Pic
* With infusion



H-2150 Subretinal Pic



H-2115 Membrane Spatula



H-2710 Scleral Depressor



H-2100 Membrane Peeler

Howard Instruments, Inc.

Website: howardinstruments.com
email: howard@howardinstruments.com
Phone: (205) 553-4453 Fax: (205) 556 9267
4749 Appletree
Tuscaloosa, AL 35405

