

## Microsurgical Accessories

A variety of products to complement ophthalmic microsurgery and minor office procedures.

### Fluid Absorption and Administration



**9001 Eye Spear**  
Fast wicking action for quick absorption



**9003 Absorbent Stick - Sterile**  
**9004 Absorbent Stick - Non-Sterile**  
2 <sup>5</sup>/<sub>8</sub> in (6.6cm) overall length <sup>3</sup>/<sub>16</sub> in (5mm) diameter, 45° angled ends  
Soft, non-woven material is rigid enough to withstand gentle pressure. Both ends may be used for absorption




**9006 Surgical Sponge - Gills**  
2 x 8mm - Lint free PVA, useful for applying or removing fluid in a precise location



**9101 Instrument Wipe**  
3 x 3 in (7.5 x 7.5cm) Lint free PVA for cleaning delicate instruments. Soft when moistened and will not drag on blades or tips




**9301 .20 Micron Filter**  
Filters out particles and microbiological contamination .22 micron and larger.  
Male/female luer lock fittings



**9302 Tubing Extension Set**  
Length: 10 in (25cm) ID: .04 in (1mm)  
OD: .08 in (2mm)  
Male/female luer lock fittings




**9303 Male to Male Connector**  
Overall length 29mm  
Luer slip fittings connect instrument to tubing  
Especially useful for bimanual I and A to connect cannula to I/A tubing



**9701 Incise Film**  
Overall size: 5 x 9 in (12.7 x 22.8cm)  
Adhesive area: 5 x 7 <sup>1</sup>/<sub>4</sub> in (12.7 x 18.4cm)  
Small incise covers lashes, lids and brow

**9802 Eye Fluid Wick**  
Highly absorbent lint free PVA wick  
Removes fluid and directs it away from surgical area



**9803 Surgical Absorbent Kit**  
1 - 9101 Instrument Wipe  
1 - 9202 Corneal Light Shield  
1 - 9802 Eye Fluid Wick  
Popular microsurgical accessories combined into one package for convenience and cost savings