



Oculoplastic

For surgical treatment, repair and irrigation of the nasolacrimal system.


Lacrimal Intubation and DCR




- 7001 Lacrimal Intubation Set**
Stainless steel probes:
23G x 4 1/2 in (.60mm x 11 cm)
Silicone tubing
Length: 14 in (35cm)
ID: .012 in (.30mm)
OD: .025 in (.60mm)



- 7002 Lacrimal Intubation Tubing**
Silicone tubing
Length: 14 in (35cm)
ID: .012 in (.30mm)
OD: .025 in (.60mm)



- 7011 DCR Set**
Stainless steel probes:
20G x 1 3/4 in (.90 x 45mm)
Silicone rod
Length: 15 3/4 in (40cm)
Diameter: .032 in (.80mm)



- 7012 DCR Set - Angled**
Angled 45°, 15mm from end
Stainless steel probes:
20G x 1 3/4 in (.90 x 45mm)
Silicone rod
Length: 15 3/4 in (40cm)
Diameter: .032 in (.80mm)

Lacrimal Irrigation



- 7021 Lacrimal Cannula**
21G x 1 1/4 in (.80 x 33mm)



- 7026 Lacrimal Cannula**
26G x 1 1/8 in (.45 x 28mm)