



Certificate:
09 105 89060

Inami & Co., Ltd.

ATTO SERIES

INAMI SLIT LAMP

SLIT LAMP ATTO L-0391 ~ L-0394



The INAMI ATTO SERIES SLIT LAMPS are not only providing practitioners with utmost satisfaction thanks to its excellent image clearness, but also great advantage of convenient carry.



Applanation Tonometer can be attached.



Continuous Zoom(10-30X)



Convenient Built-in Transformer



Suitable Carrying Case(Optional)

Specifications of L-0391

●MICROSCOPE	
Type:	Greenough Type Stereoscopic Microscope
Magnification Changer:	2 steps by changing Objective Lens
Eyepieces:	10X(16X)
Total Magnifications:	10X, 16X(16X, 25.6X)
Field of View:	20.13mm(16.10mm)
Pupillary Adjustment:	52-90mm
Focus Distance:	100 mm
Diopter Adjustment:	±6D
●SLIT ILLUMINATION	
Slit Width:	0-10mm
Aperture Diameter:	0.3, 6.5, 10,12mm
Slit Length:	0-10mm
Slit Rotation:	±90°
Filters:	Cobalt, Red-free, Built-in frost glass
Lamp:	6V/20W Halogen Lamp(L0329HI)
●BASE	
Vertical Movement:	30mm
Longitudinal Movement:	70mm
Lateral Movement:	115mm
Fine Base Movement:	±10mm
Input Voltage:	AC.100~240V 50/60Hz Built-in transformer
Power Consumption:	50VA
●CHINREST	
Vertical Movement:	75mm
●DIMENSION & WEIGHT	
Weight:	7kgs
Dimension:	240x355x480



L - 1BW
GREENOUGH TYPE MICROSCOPE



OCL -006
EYE PIECES 10X



OCL -002
EYE PIECES 16X



L-5180 Applanation Tonometer

Specifications of L-0392

●MICROSCOPE	
Type:	Galilean-Type Binocular Microscope
Magnification Changer:	3 Steps by drum rotation
Eyepieces:	15X
Total Magnifications:	10, 16, 25X
Field of View:	28.5, 14, 7mm
Pupillary Adjustment:	48-78mm
Focus Distance:	100 mm
Diopter Adjustment:	±6D



OCL -009
EYE PIECES 15X



MS-51
3 STEPS MAGNIFICATION
(10, 16, 25X)

Specifications of L-0393

●MICROSCOPE

Type:	Galilean-Type Binocular Microscope
Magnification Changer:	5 Steps by drum rotation
Eyepieces:	12.5X
Total Magnifications:	6, 10, 16, 25, 40X
Field of View:	38, 24, 15, 9, 6mm
Pupillary Adjustment:	48-78mm
Focus Distance:	100 mm
Diopter Adjustment:	+5D ~ -6D



L-5000



L-5500
Convenient Carrying Case



MS-4
5 STEPS MAGNIFICATION
(6, 10, 16, 25, 40X)



OCL-004
EYE PIECES 12.5X



Specifications of L-0394

●MICROSCOPE

Type:	Galilean-Type Binocular Microscope
Magnification Changer:	Continuous Zoom by drum rotation
Eyepieces:	15X
Total Magnifications:	10X - 30X
Field of View:	22mm-8mm
Pupillary Adjustment:	48-78mm
Focus Distance:	100 mm
Diopter Adjustment:	±6D



ZS-3-SL
ZOOM MAGNIFICATION
(10X-30X)



OCL-009
EYE PIECES 15X



●Specifications subject to change without prior notice.



Inami & Co., Ltd.

BIJYOGI RESEARCH INSTITUTE

No. 3-14-13, Bijyogi, Toda-City,
Saitama-Pref, 335-0031, Japan
Tel. 81-48-422-0173 Fax. 81-48-421-3838
Website : <http://www.inami.co.jp>
E-Mail : ibri@inami.co.jp
®TRADE MARK OF INAMI & CO., LTD.

HEAD OFFICE

No. 24-2, Hongo 3-Chome,
Bunkyo-ku, Tokyo 113-0033, Japan
Tel. 81-3-3814-5916 Fax. 81-3-5684-2126
Website : <http://www.inami.co.jp>
E-Mail : corpinfo@inami.co.jp
®TRADE MARK OF INAMI & CO., LTD.

ZCP/1000