



***Inami & Co., Ltd.***



# **PHAROS SERIES**

## **INAMI SLIT LAMP**

**SLIT LAMP PHAROS L-0185/0187/0189/0240**



***The INAMI PHAROS SERIES SLIT LAMPS provides practitioners with the utmost satisfaction thanks to its excellent image clearness and fine design.***

**L -0185 INAMI PHAROS  
GREENOUGH TYPE SLIT LAMP**



**TWO BUILT-IN DISC,  
FILTERS AND SLIT LENGTH**




**L0160H1  
HALOGEN BULB**

**L-0189 INAMI PHAROS  
5 STEPS MAGNIFICATION SLIT LAMP**



**YF-57**  
SPECIAL ADAPTOR BUILT-IN  
YELLOW FILTER



**VARIOUS KINDS OF BUILT-IN FILTERS:**  
COBALT BLUE, RED-FREE, 50%ND AND  
HEAT ABSORPTION FILTERS



**L-0187 INAMI PHAROS  
3 STEPS MAGNIFICATION SLIT LAMP**



**L-0541  
BEAM SPLITTER**

**L-0240 INAMI PHAROS  
ZOOM MAGNIFICATION SLIT LAMP**



**L-0564  
TV CAMERA ADAPTER(f=50mm)**



### Specifications of L-0185

#### ● MICROSCOPE

Type:	Greenough Type Stereoscopic Microscope
Magnification Changer:	2 steps by changing Objective Lens
Eyepieces:	10X(16X)
Total Magnifications:	10X, 16X(16X, 25.6X)
Field of View:	20.13mm(16.10mm)
Pupillary Adjustment:	52-90mm
Focus Distance:	100 mm
Diopter Adjustment:	± 6D

#### ● SLIT ILLUMINATION

Slit Width:	0-10mm
Aperture Diameter:	0.2, 1, 3, 4, 6, 10mm
Slit Length:	0-10mm
Slit Rotation:	± 90°
Filters:	Cobalt blue, Heat absorption, Red-free, 1/2ND
Slit Inclination:	5° 10° 15° 20°
Lamp:	12V/50W Halogen Bulb(L0160H1)

#### ● BASE

Vertical Movement:	30mm
Longitudinal Movement:	70mm
Lateral Movement:	115mm
Fine Base Movement:	± 10mm

#### ● CHINREST

Vertical Movement:	75mm
Fixation Target:	Tungsten lamp 6V 0.15A
Input Voltage:	AC.100, 110, 220, 240V 50/60Hz
Power Consumption:	50VA

#### ● DIMENSION & WEIGHT

Weight:	20.5kgs
Dimension:	550(W)x350(D)x500(H)



**L-1BW**  
GREENOUGH TYPE  
MICROSCOPE



**OCL-006**  
EYE PIECES 10X



**OCL-002**  
EYE PIECES 16X



### Specifications of L-0189

#### ● MICROSCOPE

Type:	Galilean-Type Binocular Microscope
Magnification Changer:	5 Steps by drum rotation
Eyepieces:	12.5X
Total Magnifications:	6, 10, 16, 25, 40X
Field of View:	38,24,15,9,6mm
Pupillary Adjustment:	48-78mm
Focus Distance:	100 mm
Diopter Adjustment:	- 6D ~ + 5D



**MS-4**  
5 STEPS  
MAGNIFICATION  
(6, 10, 16, 25, 40X)



**OCL-004**  
EYE PIECES 12.5X

## Specifications of L-0187

### ● MICROSCOPE

Type:	Galilean-Type Binocular Microscope
Magnification Changer:	3 Steps by drum rotation
Eyepieces:	15X
Total Magnifications:	10, 16, 25X
Field of View:	28.5, 14, 7mm
Pupillary Adjustment:	48-78mm
Focus Distance:	100 mm
Diopter Adjustment:	± 6D



**OCL-009**  
EYE PIECES 15X



**MS-51**  
3 STEPS MAGNIFICATION  
(10, 16, 25X)



## Specifications of L-0240

### ● MICROSCOPE

Type:	Galilean-Type Binocular Microscope
Magnification Changer:	Continuous Zoom by drum rotation
Eyepieces:	15X
Total Magnifications:	10X - 30X
Field of View:	22mm-8mm
Pupillary Adjustment:	48-78mm
Focus Distance:	100 mm
Diopter Adjustment:	± 6D



**ZS-3-SL**  
ZOOM  
MAGNIFICATION  
(10X-30X)



**OCL-009**  
EYE PIECES 15X

*"A Passion for a Nice View"*

● Specifications subject to change without prior notice.



**Inami & Co., Ltd.**

No. 24-2, Hongo 3-Chome,  
Bunkyo-ku, Tokyo 113-0033, Japan  
Tel. 81-3-3814-5916 Fax. 81-3-5684-2126  
Website : <http://www.inami.co.jp>  
E-Mail : [corpinfo@inami.co.jp](mailto:corpinfo@inami.co.jp)  
@TRADE MARK OF INAMI & CO., LTD.



Certified Management System  
DIN EN ISO 9001 : 2000  
JIS Q 9001 : 2000  
EN ISO 13485 : 2000  
ISO 13485 : 2003



ZFG1000