

# Reading magnifiers

Not only used for reading, but wherever a magnification is needed, the best known classic visual aids are reading magnifiers. Different lens constructions made from optical PMMA or glass are available in different powers and an almost unlimited variety of designs.

SEMPRAL	2.02
ROYAL SilverLine	2.03
FUNCTIONAL	2.04
BOY • SIR • BOSS /	
BOY • SIR • BOSS SilverLine	2.05
LADY	2.06
APLAMED / APLANAT	2.07
CLASSIC	2.08
FORTE / BIFOCUS / MINORA	2.09
ELEGANT / BELLA / LUXUS	2.10



## SEMPRAL

### SEMPRAL


Black colour,  
aspheric lightweight  
lens

### Cases

Leatherette,  
black colour

### Set

### Accessories



39D / 10x / Ø 35mm  
**Ref. S88390 (2)**

for S88390  
**Ref. S11941**



28D / 8x / Ø 35mm  
**Ref. S88350 (2)**

for S88350  
**Ref. S11941**



20D / 6x / Ø 55mm  
**Ref. S88550 (2)**

for S88550  
**Ref. S11942**



16D / 5x / Ø 60mm  
**Ref. S88600 (2)**

for S88600  
**Ref. S11943**



12D / 4x / Ø 70mm  
**Ref. S88700 (2)**

for S88700  
**Ref. S11944**



8D / 3x / Ø 90mm  
**Ref. S88900 (2)**

for S88900  
**Ref. S11945**



6D / 2,5x /  
□ 100x50mm  
**Ref. S88150 (2)**

for S88150  
**Ref. S11947**



**Box SEMPRAL**  
S88390, S88350,  
S88550, S88600,  
S88700, S88900,  
S88150  
**Ref. S21400**



stand for SEMPRAL  
reading magnifiers  
**Ref. S39600**



soft case for  
S88390, S88350  
**Ref. S11960**

soft case for  
S88550, S88600  
**Ref. S11962**

soft case for  
S88700  
**Ref. S11963**

(2): internal remark



## ROYAL SilverLine

### ROYAL SilverLine

Matt silver colour,  
handle in black colour,  
aspheric lightweight  
lens

Matt silver colour,  
handle in black colour,  
aspheric lightweight  
lens with anti-reflective  
hard coating

### Cases

Leather,  
silver colour

### Sets

### Accessories



39D / 10x / Ø 35mm  
**Ref. S89395**

39D / 10x / Ø 35mm  
**Ref. S89399**


for S89395, S89399  
**Ref. S11934**



28D / 8x / Ø 35mm  
**Ref. S89355**

28D / 8x / Ø 35mm  
**Ref. S89359**

for S89355, S89359  
**Ref. S11934**



20D / 6x / Ø 55mm  
**Ref. S89555**

20D / 6x / Ø 55mm  
**Ref. S89559**

for S89555, S89559  
**Ref. S11935**



16D / 5x / Ø 60mm  
**Ref. S89605**

16D / 5x / Ø 60mm  
**Ref. S89609**

for S89605, S89609  
**Ref. S11936**



12D / 4x / Ø 70mm  
**Ref. S89705**

12D / 4x / Ø 70mm  
**Ref. S89709**

for S89705, S89709  
**Ref. S11937**



8D / 3x / Ø 90mm  
**Ref. S89905**

8D / 3x / Ø 90mm  
**Ref. S89909**

for S89905, S89909  
**Ref. S11938**



6D / 2,5x /  
□ 100x50mm  
**Ref. S89155**

6D / 2,5x /  
□ 100x50mm  
**Ref. S89159**

for S89155, S89159  
**Ref. S11939**



**Box**  
**ROYAL SilverLine**  
S89395, S89355,  
S89555, S89605,  
S89705, S89905,  
S89155  
**Ref. S21505**



**Box**  
**ROYAL SilverLine,**  
**BOY SilverLine,**  
**SIR SilverLine**  
S89355, S89605,  
S89705, S89155,  
S90105, S90205  
**Ref. S21555**



**Box**  
**ROYAL SilverLine,**  
**ELEGANT, BELLA**  
S89705, S11937,  
S89905, S11938,  
S89155, S11939,  
S56750, S57750  
**Ref. S21100**

stand for ROYAL  
SilverLine reading  
magnifiers  
**Ref. S39600**



soft case for  
S89395, S89399,  
S89355, S89359  
**Ref. S11960**

soft case for  
S89555, S89559,  
S89605, S89609  
**Ref. S11962**

soft case for  
S89705, S89709  
**Ref. S11963**



## FUNCTIONAL

### FUNCTIONAL

Grey colour,  
aspheric lightweight  
lens

### Cases

Leatherette,  
black colour

### Set

### Accessories



39D / 10x / Ø 35mm  
**Ref. S83390**



for S83390  
**Ref. S11941**



**Box FUNCTIONAL**  
S83350, S83550,  
S83600, S83700,  
S83100, S83170  
**Ref. S21050**



soft case for  
S83390, S83350  
**Ref. S11960**



28D / 8x / Ø 35mm  
**Ref. S83350**

for S83350  
**Ref. S11941**


soft case for  
S83550, S83600  
**Ref. S11961**



20D / 6x / Ø 55mm  
**Ref. S83550**

for S83550  
**Ref. S11942**

soft case for  
S83700  
**Ref. S11963**



16D / 5x / Ø 60mm  
**Ref. S83600**

for S83600  
**Ref. S11943**



12D / 4x / Ø 70mm  
**Ref. S83700**

for S83700  
**Ref. S11944**



6D / 2,5x / Ø 100mm  
**Ref. S83100**

for S83100  
**Ref. S11946**



8D / 3x /  
□ 100x75mm  
**Ref. S83170**

for S83170  
**Ref. S11948**



**BOY • SIR • BOSS / BOY • SIR • BOSS SilverLine**

**Folding magnifiers, black colour**

Folding handle, folding stand with adjustable inclination, aspheric lightweight lens



BOY, 8D / 3x /  
□ 75x50mm  
**Ref. S90100**

SIR, 6D / 2,5x /  
□ 100x50mm  
**Ref. S90200**

BOSS, 8D / 3x /  
□ 100x75mm  
**Ref. S90300**

**Folding magnifiers, SilverLine**

Folding handle, folding stand with adjustable inclination, aspheric lightweight lens



BOY, 8D / 3x /  
□ 75x50mm  
**Ref. S90105**

SIR, 6D / 2,5x /  
□ 100x50mm  
**Ref. S90205**

BOSS, 8D / 3x /  
□ 100x75mm  
**Ref. S90305**

**Cases**

Leatherette, black colour



for S90100, S90105  
**Ref. S11901**

for S90200, S90205  
**Ref. S11902**

for S90300, S90305  
**Ref. S11903**

**Folding magnifiers, grey colour**

Folding handle, folding stand with adjustable inclination, aspheric lightweight lens




BOY, 8D / 3x /  
□ 75x50mm  
**Ref. S90400**

SIR, 6D / 2,5x /  
□ 100x50mm  
**Ref. S90500**

BOSS, 8D / 3x /  
□ 100x75mm  
**Ref. S90600**

**Cases**

Leatherette, grey colour



for S90400  
**Ref. S11904**

for S90500  
**Ref. S11905**

for S90600  
**Ref. S11906**

**Sets**



**Box**  
**BOY, SIR, BOSS**  
S90100, S11901, S90200, S90300, S90500, S11905, S90600, S11906  
**Ref. S21200**



**Box**  
**ROYAL SilverLine, BOY SilverLine, SIR SilverLine**  
S89355, S89605, S89705, S89155, S90105, S90205  
**Ref. S21555**



## LADY


### LADY

Transparent mount,  
handle in black colour,  
aspheric lightweight  
lens

**Cases**  
Leather,  
black colour

### Set

### Accessories



10D / 3,5x / Ø 75mm  
**Ref. S80750**

for S80750  
**Ref. S11910**



6D / 2,5x / Ø 100mm  
**Ref. S80100**

for S80100  
**Ref. S11911**



6D / 2,5x /  
□ 100x50mm  
**Ref. S80150**

for S80150  
**Ref. S11912**



5D / 2,25x /  
□ 120x60mm  
**Ref. S80160**

for S80160  
**Ref. S11913**

**Box**  
**LADY, BIFOCUS**  
S80750, S80100,  
S80150, S53750,  
S53150  
**Ref. S21030**

stand for LADY  
reading magnifiers  
**Ref. S39600**



soft case for  
S80750, S80150  
**Ref. S11963**



## APLAMED / APLANAT

### APLAMED

Transparent mount,  
handle in white colour,  
aplanatic lightweight  
lens

### Cases

Leather,  
white colour

### APLANAT

Metal mount,  
handle in black colour,  
aplanatic glass lens

### Cases

Leather,  
black colour

### Set

### Accessories



12D / 4x / Ø 65mm  
**Ref. S74650**

8D / 3x / Ø 90mm  
**Ref. S74900**



for S74650  
**Ref. S11929**

for S74900  
**Ref. S11931**



12D / 4x / Ø 65mm  
**Ref. S61650**

10D / 3,5x / Ø 75mm  
**Ref. S61750**

6D / 2,5x / Ø 90mm  
**Ref. S61900**



for S61650  
**Ref. S11924**

for S61750  
**Ref. S11925**

for S61900  
**Ref. S11926**



**Box APLAMED,  
APLANAT, CLASSIC,  
VISOLETTE,  
magnifying rulers**  
S74650, S74900,  
S61650, S51750,  
S93064, S43214,  
S43215, S43220,  
S43221  
**Ref. S21300**



stand for APLAMED  
reading magnifiers  
**Ref. S39600**



soft case for  
S74650, S61650  
**Ref. S11962**

soft case for  
S74900, S61750  
**Ref. S11963**



## CLASSIC

### CLASSIC


Metal mount,  
handle in black colour,  
biconvex glass lens

### Cases

Leather,  
black colour

### Sets

### Accessories



10D / 3,5x / Ø 65mm  
**Ref. S51600**




for S51600  
**Ref. S11917**



**Box**  
**CLASSIC, BIFOCUS**  
S51600, 2x S51750,  
S51150, S53150  
**Ref. S21020**



stand for CLASSIC  
reading magnifiers  
**Ref. S39600**



6D / 2,5x / Ø 65mm  
**Ref. S51650**

for S51650  
**Ref. S11917**



**Box APLAMED,**  
**APLANAT, CLASSIC,**  
**VISOLETTE,**  
**magnifying rulers**  
S74650, S74900,  
S61650, S51750,  
S93064, S43214,  
S43215, S43220,  
S43221  
**Ref. S21300**



soft case for  
S51840  
**Ref. S11962**



6D / 2,5x / Ø 75mm  
**Ref. S51750**

for S51750  
**Ref. S11918**



4D / 2x / Ø 90mm  
**Ref. S51900**

for S51900  
**Ref. S11919**



4D / 2x / Ø 100mm  
**Ref. S51100**

for S51100  
**Ref. S11920**



6D / 2,5x /  
□ 80x40mm  
**Ref. S51840**

for S51840  
**Ref. S11921**



4D / 2x /  
□ 100x50mm  
**Ref. S51150**

for S51150  
**Ref. S11922**



**FORTE / BIFOCUS / MINORA**

**FORTE**

Metal mount,  
handle in black colour,  
biconvex glass lens

**Cases**

Leather,  
black colour

**BIFOCUS**

Metal mount,  
handle in black colour,  
biconvex glass lens  
with bifocal lens

**MINORA**

Metal mount,  
handle in black colour,  
biconcave glass lens

**Cases**

Leather,  
black colour

**Sets and accessories**



12D / 4x / Ø 65mm  
**Ref. S55650**



for S55650  
**Ref. S11924**



6D / 2,5x / Ø 75mm,  
16D / 5x / Ø 20mm  
**Ref. S53750**



-6D / -2,5x / Ø 75mm  
**Ref. S54750**



for S53750, S54750  
**Ref. S11918**



**Box**  
**LADY, BIFOCUS**  
S80750, S80100,  
S80150, S53750,  
S53150  
**Ref. S21030**

10D / 3,5x / Ø 75mm  
**Ref. S55750**

for S55750  
**Ref. S11925**

4D / 2x /  
□ 100x50mm,  
12D / 4x / Ø 30mm  
**Ref. S53150**

for S53150  
**Ref. S11922**



**Box**  
**CLASSIC, BIFOCUS**  
S51600, 2x S51750,  
S51150, S53150  
**Ref. S21020**

4D / 2x / Ø 120mm  
**Ref. S55120**



stand for BIFOCUS  
reading magnifiers  
**Ref. S39600**



soft case for  
S55650, S55750,  
S53750, S53150,  
S54750  
**Ref. S11963**



**ELEGANT / BELLA / LUXUS**

**ELEGANT**

Silver mount,  
wooden handle,  
biconvex glass lens

**BELLA**

Golden mount,  
wooden handle,  
biconvex glass lens

**LUXUS**

Golden mount,  
handle with  
leather trimming,  
biconvex glass lens

**Cases**

Leather,  
brown colour

**Set**

**Accessory**



6D / 2,5x / Ø 75mm  
**Ref. S56750**



6D / 2,5x / Ø 75mm  
**Ref. S57750**



6D / 2,5x / Ø 75mm  
**Ref. S58750**



for S56750, S57750,  
S58750  
**Ref. S11932**



**Box**  
**ROYAL SilverLine,**  
**ELEGANT, BELLA**  
S89705, S11937,  
S89905, S11938,  
S89155, S11939,  
S56750, S57750  
**Ref. S21100**



soft case for  
S56750, S57750,  
S58750, S58150  
**Ref. S11963**

4D / 2x /  
□ 100x50mm  
**Ref. S58150**

for S58150  
**Ref. S11933**