

Maloney Hand-Held Surgical Keratoscope

Coaxial light from the microscope is directed through the center of the instrument to display rings on the cornea.

28-4430



Van Loenen Cylindrical Keratoscope

Side illuminating the barrel of the instrument displays rings on the cornea.

JEDMED/Van Loenen Keratoscope, non-sterile, pk/5 28-4420


JEDMED/Van Loenen Keratoscope, non-sterile, each 28-4421



The **McIntyre Oculo-Pressor™**

A high density alloy, compression molded for compactness, provides constant and gentle force equivalent to 30mm of Hg pressure.

29-0450



Honan Intraocular Pressure Reducer (Honan's Balloon)

Permits surgeons to effectively lower intraocular pressure by a known low external pressure prior to intraocular surgery.

29-0500

