

# *TL-90M / 134M*

*Trial lens sets*



**MAGNON**

# TL-90M / 134M

**MAGNON**

## Trial lens sets

The MAGNON TL-90M and TL-134M are the economical trial lens sets with metal rims to meet minimum requirement.

### FEATURES

- ★ Double curve lenses with 38mm in diameter
- ★ Metal rims (Gold-plated for cylinder lenses)
- ★ Space for MAGNON trial frame
- ★ Portable leatherette carrying case

### SPECIFICATIONS

	TL-90M	TL-134M
Spheres	17 pairs each of plus and minus	24 pairs each of plus and minus
Cylinders	8 pairs of minus only	12 pairs of minus only
Prisms	4 pieces	8 pieces
Accessories	2 pieces	6 pieces
	● : Occluder .....1	● : Occluder .....1
	F : Frosted glass .....1	○ : Pin hole (1 mm) .....1
		C : Cross cylinder $\pm 0.50D$ .....1
		F : Frosted glass .....1
		N : 0.00D lens .....1
		M : Maddox .....1

#### TL-90M

- Spheres	- Cylinders	+ Spheres
10.00	3.00	10.00
8.00	2.50	8.00
6.00	2.00	6.00
5.00	1.50	5.00
4.50	1.00	4.50
4.00	0.75	4.00
3.50	0.50	3.50
3.00	0.25	3.00
2.50		2.50
2.00		2.00
1.75		1.75
1.50	Prisms	1.50
1.25	1° 2° 3° 4°	1.25
1.00		1.00
0.75	Accessories	0.75
0.50		0.50
0.25	● F	0.25

#### TL-134M

- Spheres	- Cylinders	+ Spheres
12.00	4.00	12.00
10.00	3.50	10.00
9.00	3.00	9.00
8.00	2.50	8.00
7.00	2.00	7.00
6.00	1.75	6.00
5.00	1.50	5.00
4.50	1.25	4.50
4.00	1.00	4.00
3.75	0.75	3.75
3.50	0.50	3.50
3.25	0.25	3.25
3.00		3.00
2.75		2.75
2.50	Prisms	2.50
2.25		2.25
2.00	4° 5°	2.00
1.75	2° 3°	1.75
1.50	1°	1.50
1.25	0.5°	1.25
1.00		1.00
0.75	Accessories	0.75
0.50		0.50
0.25	● ○ CFNM	0.25

Specifications and design are subject to change for improvement without notice.

**H. OGINO & CO., LTD.**

Yokohama Port P.O.Box 170, Yokohama 221-8691, Japan  
Tel: (045)862-5261 Fax: (045)851-5467

020102A Printed in Japan