

# *TLX-232*

*Trial lens set*



**MAGNON**

# TLX-232

## Trial lens set

# MAGNON

The MAGNON TLX-232 is corrected curve meniscus trial lens set covering most of lens powers which doctors require. The corrected curve additive design compensates for thickness and space between combined lenses, which allows the exact correction to be read directly from the total lens power in the trial frame.

### FEATURES

- ★ Meniscus lenses of 25mm in diameter
- ★ Color coded and anodized aluminum rims of 38mm in diameter
- ★ One side marked lens power for insuring correct position in the trial frame
- ★ Light weight carrying case
- ★ Plastic tray with space for the trial frame

### SPECIFICATIONS

Spheres	35 pairs each of plus and minus		
Cylinders	18 pairs each of plus and minus		
Prisms	10 pieces		
Accessories	10 pieces		
	B.L : Occluder .....1	R	: Red lens .....1
	O.O : O.O lens .....1	G	: Green lens .....1
	P.H : Pin hole (1mm) .....1	R.M	: Red maddox .....1
	SS : Silt (1mm) .....1	F.G	: Frosted glass .....1
	CL : Test mark for PD .....2		

- Spheres	Prisms		Accessories		+ Spheres
20.00	10	12	R.M	F.G	20.00
18.00	6.0	8.0	R	G	18.00
16.00	4.0	5.0	CL	CL	16.00
15.00	2.0	3.0	P.H	S.S	15.00
14.00	0.5	1.0	B.L	0.00	14.00
13.00					13.00
12.00					12.00
11.00					11.00
10.00					10.00
9.00					9.00
8.00					8.00
7.50					7.50
7.00					7.00
6.50					6.50
6.00					6.00
5.50					5.50
5.00					5.00
4.50					4.50
4.00					4.00
3.75					3.75
3.50					3.50
3.25					3.25
3.00					3.00
2.75					2.75
2.50					2.50
2.25					2.25
2.00					2.00
1.75					1.75
1.50					1.50
1.25					1.25
1.00					1.00
0.75					0.75
0.50					0.50
0.25					0.25
0.12					0.12

Specifications and design are subject to change for improvement without notice.

### H. OGINO & CO., LTD.

Yokohama Port P.O.Box 170, Yokohama 231-0891, Japan  
Tel (045)662-5261 Fax (045)651-5467