



Novamed Ceramisys



Novamed Ceramisys Ltd brings you the very latest in
Bioceramic and Hydroxyapatite Orbital Implant technology

Distributors required in some countries.

***Until 31st March 2007 buy 5 implants and receive
an additional one free of charge!***


- Pure Alumina, not recognised as a foreign body and does not illicit immune response.
- No possibility of resorption and no ion release
- The implant is supplied sterile in a range of sizes 14, 16, 18, 19, 20, 21 and 22mm
- Ultra-porosity, typically 200-800um pores
- 100% interconnected with no dead pockets
- Ultra lightweight, over 30% lighter than reported for alternative products
- We also supply Hydroxyapatite Orbital Implants in a similar range of sizes
- **Very competitive pricing**
Less expensive than other similar products

***For further information or details of
Distribution in your country please call or
email sales@novamed.co.uk***

CE 95132



*Novamed Ceramisys
Alumina orbital implant*



*Human cancellous Bone
(from Tibial Plateaux)*

Novamed Ceramisys Ltd
PO Box 232, Montague Street
Dundee, DD5 2YN, Scotland
Tel : + 44 (0) 1382 774777
Fax : + 44 (0) 1382 775777
Email : sales@novamed.co.uk
Web: www.novamed.co.uk