

Binocular Indirect Ophthalmoscopy (BIO) Lenses



14 Diopter Lens

High magnification for detailed examination of macula and optic disc.

Product Code	Field of View	Magnification
01-14	37°	4.29x



20 Diopter Lens

Most common lens for B.I.O.
High resolution image.

Product Code	Field of View	Magnification
01-20	50°	2.97x

Indirect Lenses



Triple Two Panfundus Lens

22D lens for general fundus exam with the binocular indirect ophthalmoscope. Large diameter and unique optical design combine magnification with very wide field.

Product Code	Field of View	Magnification
01-222	60°	2.72x



28 Diopter Lens

Excellent general purpose lens. Small diameter easy to handle. Popular for examining children.

Product Code	Field of View	Magnification
01-28	58°	2.13x