

Securefit™ flange for ultimate control.



Fundus Diagnostic Lens

This "Goldmann" type fundus lens provides a direct image of the posterior pole.



Product Code	Contact Diameter	Lens Height	Field of View	Lens Magnification
OGF	15mm	16.5mm	36°	.93x
OGF-2*	15mm	16.5mm	36°	.93x



Magna View Gonio Lens

The best lens available for gonioscopy and laser trabeculoplasty. One 62° mirror. Tilted anterior surface corrects image and laser beam astigmatism. Unsurpassed resolution. Laser beam spot size magnification .77x.



Product Code	Contact Diameter	Lens Height	Lens Magnification
OMVGL	15mm	23.5mm	1.3x