



Aladdin Classic

- > integrated, one-piece design
- > 14" screen for easy reading and writing
- > adjustable contrast reduces glare to sensitive eyes
- > American Society of Aging award winning design

Remember reading for pleasure? With the Aladdin Classic, you can return to activities you enjoy most—reading books and writing letters—even doing close handwork. The oversized reading tray accommodates newspapers, magazines, even large photo albums. Plus, with magnification to 50X, you'll be able to attend to personal business—like viewing medicine bottles and balancing your checkbook.

A practical all-in-one design and easy-to-use controls are just two features that earned the Classic a Gold Award from the American Society on Aging.

Different viewing modes and features help you enjoy a range of activities:

- > For longer reading sessions, white letters on a low-glare black background.
- > Photo viewing mode puts snapshots, illustrations and 3D objects into clear focus.
- > Sliding picture control settings allow you to adjust to new tasks quickly.
- > Margin stops make it easier to read columns of text.

Telesensory > is the leading provider of magnification systems for people with low vision. We have restored independent living to more than 200,000 satisfied customers.



For more information or a free no-obligation demonstration, call **1-800-804-8004**.

Features > Aladdin Classic

Magnification from 4.5X to 50X. Newsprint can be magnified to over 3" high.

No assembly required / Single piece construction

Easy-to-Use Controls:

- > Size adjustment lever
- > Sliding picture control adjustment
- > Focus control knob

Display > Black & White 14" monitor, 12.8" viewing area

Multiple Viewing Options for Different Needs:

Bright Negative Viewing Mode > uses bright white letters on a low-glare black background for maximum reading comfort.

Bright Positive Viewing Mode

Photo Viewing Mode > is a medium contrast setting ideal for viewing photos, illustrations or 3-dimensional objects, such as a prescription bottle, in fine detail.

Reading Table > a large non-glare, platform for reading oversized books that moves freely in any direction—or locks into place. Adjustable margin stops make it easier to read text in columns.

Money-back Guarantee > Telesensory products are backed by a 30 day money-back guarantee (less shipping) and a two year manufacturer's limited warranty against defects in material and workmanship.

Information and Sales > Representatives are available Monday through Friday from 7 am to 4 pm PST to answer questions and take your orders.

Specifications

Magnification: 4.5X to 50X

Unit Dimensions:

Width 15.6" (39.6 cm)

Depth 21.5" (54.6 cm)

Height 22" (55.9 cm)

Monitor: Black & White, 14" (35.6 cm), 12.8" (32.5 cm) viewing area (diagonal)

Weight: 37 lbs (16.8 kg)

Power Requirements:

AL2A: 60 HZ, 115 VAC ($\pm 10\%$), 50W max. AL5A: 47-63 HZ, 100-240 VAC ($\pm 10\%$) 50W max. (Unit includes removable power cord.)

Fuse Types: 2 amp time delay, 250V

Environment:

Operating Temperature: +41F to +122F (+5C to +50C)

Storage Temperature: -4F to +140F (-20C to +60C)

Approvals:

AL2A: c-ETL-us and FDA,
AL5A: c-ETL-us, CE, GS, VCCI, SEMKO, DEMKO, NEMKO, FIMKO and FDA,
AL5A-MPRII: GS, CE, SEMKO, DEMKO, NEMKO, FIMKO, MPRII and FDA

Accessories

Carrying Case: This sturdy case protects your Classic, and allows for easy transport.

Dust Cover

Contact InSiPhil at :

Toll free in the US: 1800 - 804 - 8004 or 408 - 616 - 8700

Address : Plo 160, Jalan Cyber 2, Senai Industrial Estate III
81400 Senai, Johor, Malaysia. Tel : (60)7 598 5010

www.insiphil.com / www.telesensory.com

email : info@insiphil.com