



## Aladdin Apex

- > adjustable 17" or 20" flat panel display
- > integrated, one-piece design
- > auto focus with manual override
- > full color with 24 color select options

**NEW!**  
Flat Panel



Part of the Aladdin one-piece desktop family, the Aladdin Apex incorporates an adjustable flat panel monitor with the same color and professional features you've come to rely on.

Enjoy reading books and magazines, doing crossword puzzles and examining hard-to-see objects such as pill bottles — all with the ease of the Aladdin Apex's adjustable flat panel screen. With the Apex you can raise or lower the screen, bringing the image to your eye level. Or slide the screen forward and tilt to enjoy a glare-free image at your optimum reading distance.

Reading is made easy with this full-featured, integrated flat panel system:

- > 24 color select options offer high contrast
- > Keep your place with line markers and shadow masking
- > Continuous auto focus with manual override
- > Monitor folds down for easy storage

**Telesensory** > With over 35 years of success, Telesensory continues to be the company people trust and depend on for the highest quality low vision solutions and dedicated personal customer service.



**Telesensory**  
*focused on you™*

For more information or a free no-obligation demonstration, call **1-800-804-8004**.

## Features > Aladdin Apex

Magnification from 3X to 48X for the 17" model; 3.5X to 56X for the 20" model. Newsprint can be magnified to over 4" high.

No assembly or setup required.

Easy-to-Use Controls:

- > Size adjustment lever
- > Picture control with three zones
- > Continuous auto-focus with manual override
- > Shadow mask/Line markers
- > Variable brightness

Display > Full color 17" or 20" flat panel monitor

Multiple Viewing Options for Different Needs:

Adjustable Flat Panel Display > lowers and extends to accommodate each person's distinct height and eye level preferences. Tilting the display helps reduce glare.

Color Select Mode > lets you view reading materials in 24 different color combinations to accommodate specific vision needs.

Negative Text Viewing Mode > uses bright white letters on a low-glare black background for maximum reading comfort.

Two Full Color Viewing Modes > reading mode offers high-quality color contrast for easy viewing of materials with photos and text. Photo mode is a medium contrast setting ideal for viewing photos, illustrations or 3-dimensional objects in detail.

Reading Table > a non-glare platform for reading that moves freely in any direction—or slowly with a table brake. Adjustable margin stops make it easier to read text in columns.

Money-back Guarantee > Telesensory products are backed by a 30 day money-back guarantee (less shipping) and a two year manufacturer's limited warranty against defects in material and workmanship.

Information and Sales > Representatives are available Monday through Friday from 7 am to 4 pm PST to answer questions and take your orders.

## Contact InSiPhil at :

Toll free in the US: 1800 - 804 - 8004 or 408 - 616 - 8700  
Address : Plo 160, Jalan Cyber 2, Senai Industrial Estate III  
81400 Senai, Johor, Malaysia. Tel : (60)7 598 5010  
www.insiphil.com / www.telesensory.com  
email : info@insiphil.com

## Specifications

Magnification:

17" Model – 3X to 48X;  
20" Model – 3.5X to 56X

Unit Dimensions (folded)

17" Model:

Width: 16" (40.6 cm)

Depth: 20.5" (52 cm)

Height: 18.5" (47 cm)

20" Model:

Width: 20" (50.8 cm)

Depth: 21.5" (54.6 cm)

Height: 18.5" (47 cm)

Monitor: Full Color, 17" (43.2 cm);  
20" (50.8 cm)

Weight:

17" Model: 47 lbs (21.4 kg)

20" Model: 52 lbs (23.5 kg)

Power Requirements:

47-63 HZ, 100-240 VAC (±10%)

85W max. (Unit includes removable power cord.)

Fuse Types: 2 amp time delay, 250V

Environment:

Operating Temperature: +41F to +122F  
(+5C to +50C)

Storage Temperature: -4F to +140F  
(-20C to +60C)

Approvals: c-ETL-us, CE, VCCI, CB  
and FCC