



## Olympia

- > portable video magnifier
- > eight magnification settings – from 4.3X to 26X
- > connects to TV monitor for even larger magnification
- > adjustable brightness



This powerful portable magnifier gives you the power of a desktop video magnifier, with the freedom of a portable product. Olympia provides maximum flexibility and control—including a wide range of magnification settings, adjustable brightness, and a tilted screen to reduce glare. Enjoy hours of reading for work, school or pleasure on the built-in 7" display. You'll want to keep the versatile Olympia with you at home or when traveling. Experience even greater magnification levels by connecting Olympia to any compatible TV monitor. And, we've included a writing stand so you can easily fill out forms, complete your homework, and write letters and checks.

Innovative viewing features and flexible settings:

- > Our wide screen view and tilted screen make it easier to read and write in a variety of circumstances.
- > Choose color, positive or negative text, or yellow on blue viewing mode and adjust brightness to meet your visual needs.

There's no assembly or set up required. With the carrying case, battery pack and optional car charger you can take it anywhere.

**Telesensory** > With over 35 years of success, Telesensory continues to be the company people trust and depend on for the highest quality low vision solutions and dedicated personal customer service.



For more information or a free no-obligation demonstration, call 1-800-804-8004.

## Features > Olympia

Magnification of 4.3X to 26X

26X magnification increases newsprint to 2.0" high.

Magnification range from 12X to 76X with 20" monitor.\*

No assembly or set up required.

Easy-to-Use Controls:

- > Ergonomic magnification control slide
- > Tilted screen with wide screen view
- > Positive, negative, color and yellow on blue display options
- > Adjustable brightness

Display > 7 inch (17.18 cm) diagonal, color LCD.

Tilted screen with wide screen view.

Multiple Options for Viewing and Writing:

**Four Viewing Modes >** Color Mode for viewing photos and magazines. Positive Text Mode, using black letters on a bright white background. Negative Text Mode, using white letters on a black background. Yellow letters on a blue background provide a more comfortable contrast for some eye conditions.

**Writing Feature >** Two magnification levels 4.5X and 9X with 4.75" clearance under the writing stand.

**Money-back Guarantee >** The Olympia is backed by a 30 day money-back guarantee (less shipping) and a two year manufacturer's limited warranty against defects in material and workmanship. The battery is backed by a 1-year warranty.

**Information and Sales >** Representatives are available Monday through Friday from 7 am to 4 pm PST to answer questions and take your orders.

\* Requires phono-type jack. Non-NTSC TV/monitor requires converter.

**Contact InSiPhil at :**

Toll free : 1800 - 804 - 8004 or 408 - 616 - 8700

Address : Plo 160, Jalan Cyber 2, Senai Industrial Estate III  
81400 Senai, Johor, Malaysia. Tel : (60)7 598 5010

www.insiphil.com / www.telesensory.com

email : info@insiphil.com

## Specifications

**Magnification (3 viewing options):**

Standard: 4.3X, 7.8X, 13X

Max Magnification: 8.6X, 15.6X, 26X

Writing: 4.5X, 9X

**Unit Dimensions:**

Length: 8.2" (20.8 cm)

Width: 5.4" (13.7 cm)

Height: 4.25" (10.8 cm)

**Weight (Excluding Battery):**

2.5 lbs (1.14 kg)

**Battery:**

Use time: 2.5 to 6 hours

Charge time: 2 to 5 hours

**Display:** Color, 7" (17.8 cm)

diagonally (3.5" H x 6.1" W)

(8.9 cm H x 15.6 cm W)

**Power Requirements:** 47-63 HZ,  
90-264 VAC, 3.4A max load or 12V DC,  
1.5A to 3.5A charging battery pack.

**Environment:**

Operating temperature: +41F to +122F  
(+5C to +50C)

Storage temperature: -4F to +140F  
(-20C to +60C)

## Accessories

**Carrying Case (included):** This sturdy black case protects the Olympia and battery pack and makes travel easy.

**Battery Pack (included) and Car Adapter (optional):** The Olympia external battery pack provides portable power. Recharge the battery in just two hours at home—or on the road with our handy car adapter.

**Writing Stand (included):** Specially designed with space for your hand and maximum brightness, this folding stand provides a stable platform to write letters, checks and do paperwork.