



Pico

NEW!

- > pocket-sized video magnifier
- > integrated rechargeable battery
- > four viewing modes
- > magnification from 3X to 11X



Read the fine print anytime, anywhere! Pico fits easily in your purse or pocket, and the integrated battery makes it easy to carry wherever you go. With Pico, you can even read in dimly-lit areas. Pico offers yellow on blue and positive text modes to accommodate specific vision needs. The bright display and viewing angle provide precise viewing and comfortable positioning.

Created for convenience, Pico adds to Telesensory's line of portable video magnifiers. Pico is so small you will always want it with you—shopping for greeting cards, looking at price tags in stores and reading food ingredients in the supermarket. Easily view three-dimensional objects such as pill bottles on Pico's built-in color display. Fold-down legs allow you to vary magnification, write and sign checks and tilt the screen toward you for reading.

Convenience is at your fingertips:

- > Read bills, schedules, grocery lists, receipts, menus and more in bright color.
- > Take it to the bank to read ATM instructions and transaction records.
- > Sign checks and bills, or write a quick note even when you're on the go.

Living with low vision just got easier. Wherever life takes you, equip yourself with a personal reader that goes the distance.

Telesensory > With over 35 years of success, Telesensory continues to be the company people trust and depend on for the highest quality low vision solutions and dedicated personal customer service.



For more information or to purchase a product, call **1-800-804-8004**.

Features > Pico

Magnifies print and photos from 3X to 11X.

5X magnification increases newsprint to 0.4" high.

No assembly or set up required.

Easy-to-Use Controls:

- > One button to turn Pico on and off
- > One button to select color, positive or negative text, yellow on blue viewing modes
- > Lighting control switch to dim or turn off LED lights

Display > 4 inch (10.16 cm) diagonal, color LCD

Multiple Options for Viewing and Writing:

Four Viewing Modes > Color Mode for viewing photos and magazines. Negative text mode provides high contrast white letters on a black background. Positive text mode provides black letters on a bright white background. Yellow letters on a blue background provide a more comfortable contrast for some eye conditions.

Versatile Fold-down Legs > The versatile folding legs provide different magnification levels. Both legs unfolded—5X magnification; top leg unfolded—11X; bottom leg unfolded—8X. Lift Pico off surface to reduce magnification to 3X. Plus, fold-down legs allow you to write, sign checks or credit receipts.

Rechargeable Battery > Plug the AC Adapter into Pico to re-charge the unit. You can even read while you charge without increasing the time it takes to charge the unit.

Money-back Guarantee > The Pico is backed by a 30-day money-back guarantee (less shipping) and a two year manufacturer's limited warranty against defects in material and workmanship. The battery is backed by a two year warranty.

Information and Sales > Representatives are available Monday through Friday from 7 am to 4 pm PST to answer questions and take your orders.

Contact InSiPhil at :

Toll free in the US: 1800 - 804 - 8004 or 408 - 616 - 8700
Address : Plo 160, Jalan Cyber 2, Senai Industrial Estate III
81400 Senai, Johor, Malaysia. Tel : (60)7 598 5010
www.insiphil.com / www.telesensory.com
email : info@insiphil.com

Specifications

Magnification: 3X to 11X

Unit Dimensions

Length: 6.5" (16.51 cm)

Width: 3.5" (8.89 cm)

Height: 1" (2.54 cm)

Weight: 10 ounces (0.28 kg)

Battery:

Use time: approx. 1 hour, 15 minutes

Charge time: approx. 2 hours

Display: Color, 4" (10.6 cm) diagonally (2.4" H x 3.2" W) (6.1 cm H x 8.1 cm W)

Power Requirements:

Input: 50-60 HZ, 100-240 VAC

Battery Charging: 10V DC; 0.8A (max)

Battery Charging and Reading: 10V DC; 1.4A (max)

Accessories

Carrying Case (included):

This custom case protects Pico and makes it easy to carry in your purse or on your belt.

Wrist Strap (included):

Designed to keep the Pico safely on your arm.

Battery Charger (included):

Plug the AC Adapter into Pico for approximately 2 hours to recharge the unit.