

# Readit

The Readit from VisionAid International Ltd is the World's most powerful desktop reading solution. The Readit will capture a page of printed text up to A3 in size and scroll the text on the screen whilst reading it aloud. The Readit also has a dedicated local and distance viewing camera meaning you can read handwriting and distant objects such as white boards. It's also a full PC!



## Readit benefits:

- Capture and recognise most documents in 30 seconds. The text can then be displayed in the following enhanced views, greatly increasing reading speed:



**Overlay View**



**Horizontal View**



**Vertical View**



**Column View**

- Scroll the recognised text with up to 60X magnification. You can also choose to view the text in 17 different colour combinations.
- Four high quality voices allow you to listen to your captured documents as well as see them enlarged on screen.
- The Readit's Live camera enables you to view handwriting (that can not be captured) and view distant objects. The live camera has 16 enhanced colour combinations as well and can magnify up to 40X.
- Readit is also a full Windows XP Home PC and you can use any of your screen reading or magnification software with it. It comprises a 19, 22 or 24" screen, 2.66 Ghz Core2Duo processor, 2gb of RAM, 250GB HDD and DVD RW.

