



High Technology for  
Your Ophthalmology



**CATALOGUE 2007**

**WWW.CILITA.COM**

**Dear Sir / Madam,**

**We are proud to present you CILITA LTD, a reliable manufacturer and an established retailer of ophthalmic surgical instruments based in Russia.**

**Possessing a big experience in manufacturing and marketing of ophthalmic products we tend to produce high quality microsurgical instrumentation for application in the field of General Ophthalmology, Cataract and Vitreoretinal Surgery.**

**Ophthalmic instrument market has been always divided between the suppliers from Europe and America which prices were always groundlessly high and from India and China which quality continues to be too poor for the growing requirements of ophthalmic practices. We have established our company to manufacture microsurgical instruments of the highest quality but to market them at moderate, cost-effective prices. Please, compare them and you will see that we offer the best price/quality ratio available in the market to meet or exceed your requirements.**

**Take a few minutes and visit our web site [www.cilita.com](http://www.cilita.com) and you will see that we offer a wide variety of ophthalmic surgical instruments.**

**CILITA LTD sells its products internationally directly through the online-shop and through independent distributors all over the world.**

**We appreciate different ways of cooperation and are looking forward to working with you.**

# BLADE HOLDERS, CALIPERS & MARKERS


## Blade Holder



## B-010 Short Model

- Round knurled handle with end lock
- Polished finish
- Overall length 120 mm

## Castroviejo Caliper



## C-010T titanium Straight

- Measures from 0 to 20 mm
- Scale reads on both sides
- Polished finish
- Overall length 87 mm

## Osher-Neumann Low-Profile Corneal Marker



## C-0304 8 lines

- Round serrated handle
- Overall length 130 mm

## Corneal Transplant Marker



## C-140

- Internal diameter 7.0 mm
- External diameter 11 mm
- Round serrated handle
- Overall length 130 mm

## LASIK Marker




## C-176

- Optical zone 8 mm
- 3 marking lines, asymmetrical, to ensure correct placement of corneal flap
- Smooth edges
- Designed for both nasal and superior hinge
- Overall length 130 mm

# ANTERIOR CHAMBER FORCEPS

Anterior Chamber Forceps, Angled



F-0290S stainless steel

- Angled shaft
- 1x2 teeth
- Flat serrated handle
- Polished finish
- Overall length 90 mm

# CAPSULORHEXIS FORCEPS

Utrata Capsulorhexis Forceps Straight Jaws, Small Size



F-0331T titanium

*Cystotome tips*

Jaws 11 mm

- Needle holder handle
- Overall length 115 mm

Utrata Capsulorhexis Forceps Straight Jaws Flat Handle



F-0301T titanium

*Cystotome tips*

Jaws 11 mm for Extracapsular  
Clear Cornea & Tunnel

- Flat serrated handle
- Overall length 110 mm

Utrata Capsulorhexis Forceps Straight Jaws Round Handle



F-0311T titanium

*Cystotome tips*

Jaws 11 mm for Extracapsular  
Clear Cornea & Tunnel

- Round serrated handle
- Overall length 107 mm

Utrata Capsulorhexis Forceps Straight Jaws Flat Handle



F-0321T titanium

*Cystotome tips*

Jaws 11 mm for Extracapsular  
Clear Cornea & Tunnel

- Flat serrated handle
- Overall length 110 mm

# CAPSULORHEXIS FORCEPS

**Utrata Capsulorrhexis Forceps, Flat Handle**



**F-0300T titanium**

***Cystotome tips***

- Jaws 11 mm for Extracapsular Clear Cornea & Tunnel
- Flat serrated handle
- Overall length 82 mm

**Side Port Capsulorhexis Forcerps 20G, Intraocular Handle**



**F-037/V-003**

- Intraocular titanium handle
- Cystotome tips
- 20 Gauge forceps

**Side Port Capsulorhexis Forcerps 23G, Intraocular Handle**



**F-0374/V-003**

***for Cold Phaco***

- Intraocular titanium handle
- Cystotome tips
- 23 Gauge forceps

**Inamura Capsulorhexis Forceps**



**F-0391**

***Cystotome tips, stainless steel***

- Jaws 10 mm
- To maintain alignment of tips, prevent leakage of viscoelastic from anterior chamber
- Tips do not open more than 2 mm
- Needle holder handle
- Overall length 115 mm

**Inamura Capsulorhexis Forceps**



**F-0392**

***Cystotome tips, stainless steel***

- Jaws 12 mm
- To maintain alignment of tips, prevent leakage of viscoelastic from anterior chamber
- Tips do not open more than 2mm
- Needle holder handle
- Overall length 115 mm

# LID & CHALAZION FORCEPS

**Snellen Entropium Forceps, Small Left**



- F-1904S stainless steel**
- Left small
  - Solid lower plate , open upper plate
  - Locking thumb screw mechanism
  - Polished finish
  - Overall length 90 mm

**Snellen Entropium Forceps, Small Right**




- F-1905S stainless steel**
- Right small
  - Solid lower plate , open upper plate
  - Locking thumb screw mechanism
  - Polished finish
  - Overall length 90 mm

**Desmarres Chalazion Forceps Oval**



- F-1906T titanium,**  
large size, 27mm x 17mm
- 4-1907T titanium,**  
medium size, 12mm x 10mm
- Solid lower plate, open upper plate
  - Locking thumb screw mechanism
  - Polished finish
  - Overall length 90 mm

**Lambert Chalazion Forceps, Small**



- F-1908T titanium**
- Round solid lower plate, open upper plate
  - 8 mm inside diameter of upper plate
  - Locking thumb screw mechanism
  - Polished finish
  - Overall length 92 mm

**Lambert Chalazion Forceps, Medium**



- F-1909T titanium**
- Round solid lower plate, open upper plate
  - 12 mm inside diameter of upper plate
  - Locking thumb screw mechanism
  - Polished finish
  - Overall length 97 mm

**Desmarres Chalazion Forceps**



- F-1912T titanium,**  
large size
- Oval solid 24 mm wide lower plate, open upper plate
  - 12mm x 20mm inside dimensions of upper plate
  - Locking thumb screw mechanism
  - Polished finish
  - Overall length 92 mm

## LID & CHALAZION FORCEPS

### Blepharoplasty Lid Clamp Forceps



**F-1920 Large size**

30 mm x 24 mm

**F-1921 Small size**

30 mm x 20 mm

- Lid clamp forceps designed for blefaroplasty to protect the globe and provide a platform for lid resection
- Locking thumb screw mechanism

## CILIA FORCEPS

### Ziegler Cilia Forceps



**F-040T titanium**

- Smooth jaws
- Serrated handle
- Polished finish
- Overall length 94 mm

### Cilia Forceps, Wide



**F-043T titanium**

- Smooth jaws
- Serrated handle
- Polished finish
- Overall length 94 mm

## CONJUNCTIVA FORCEPS

### Fechtner Conjunctiva Forceps



**F-2301T titanium**

- Delicate ring tip jaws
- Flat serrated handle
- Polished finish
- Overall length 100 mm

# CORNEAL FORCEPS

## Colibri Corneal Forceps



**F-050T titanium**  
**0.12 mm teeth**

- Flat serrated handle
- 0.12 mm, 1x2 teeth
- Tying platform
- Dull finish
- Overall length 82 mm

## Colibri Corneal Forceps



**F-0501T titanium**  
**0.12 mm teeth**

- Flat serrated handle
- 1x2 teeth
- Tying platform
- Dull finish
- Overall length 85 mm

## Hoskin Type Colibri Forceps, Pierse Tips



**F-0502T titanium**

- Flat serrated handle
- Tying platform
- Dull finish
- Overall length 82 mm
- Available in Titanium or Stainless Steel

## Colibri-Bonn Corneal Forceps



**F-0504T titanium**  
**0.12 mm teeth**

- Flat serrated handle
- 0.12 mm, 1x2 teeth
- Tying platform
- Dull finish
- Overall length 85 mm

## Micro Colibri Corneal Forceps



**F-0505T titanium**  
**0.12 mm teeth**

- Flat serrated handle
- 0.12 mm, 1x2 teeth
- Tying platform
- Dull finish
- Overall length 72 mm

## Colibri Corneal Forceps



**F-053T titanium**  
**0.12 mm teeth**

- Flat serrated handle
- 0.12 mm, 1x2 teeth
- Tying platform
- Dull finish
- Overall length 115 mm

# CORNEAL FORCEPS

## Colibri Corneal Forceps, Round Handle



**F-054T titanium**

**0.12 mm teeth**

- Round serrated handle
- 1x2 teeth
- Tying platform
- Dull finish
- Overall length 115 mm

## Straight Corneal Forceps, Bonn-Catalano Type



**F-0551T titanium**

**0.12 mm teeth**

- Round serrated handle
- 0.12 mm, 1x2 teeth
- Tying platform
- Dull finish
- Overall length 105 mm

## Catalano Corneal Forceps



**F-055T**

**1x2 teeth 0.12mm**

**F-057T**

**1x2 teeth 0.5 mm**

- Round serrated handle
- Tying platform
- Overall length 105 mm

## Bonn Corneal Forceps



**F-058T titanium**

**Medium size, 94mm**

**F-059T titanium**

**Small size, 72mm**

**0.12 mm teeth**

- Flat serrated handle
- 0.12 mm, 1x2 teeth
- Tying platform
- Dull finish

# DRESSING FORCEPS

## Straight Dressing Forceps With Delicate Serrations



### F-070T titanium

- Straight 6 mm serrated tips
- Polished finish
- Round serrated handle
- Overall length 100 mm

## Straight Dressing Forceps With Serrations



### F-071T titanium

- Straight 12 mm serrated tips
- Flat serrated handle
- Polished finish
- Overall length 100 mm

## Curved Dressing Forceps With Serrations



### F-072T titanium

- Curved serrated tips
- Flat serrated handle
- Polished finish
- Overall length 100 mm

## Adson Dressing Forceps




### F-0740T titanium

- Ideal for plastic & reconstructive surgery*
- Flat serrated handle
  - Polished finish
  - Overall length 125 mm

# FIXATION FORCEPS

## Adson Fixation Forceps



### F-0741T titanium

*Ideal for plastic & reconstructive surgery*

- 1x2 teeth 0.5 mm
- Flat serrated handle
- Polished finish
- Overall length 125 mm

## Kremer Fixation Forceps



### F-0801S stainless steel *With Lock*

- Two point fixation
- Tips separated by 13 mm
- 0.12 mm, 1x2 teeth
- Flat serrated handle
- Polished finish
- Overall length 105 mm


## Polack Corneal Transplant Fixation Forceps



### F-0814T titanium *Curved*

- Two point fixation
- Tips separated by 1.5 mm
- Pierse type tips
- Round serrated handle
- Polished finish
- Overall length 100 mm

## Castroviejo Fixation Forceps, W/O Platform



### F-0822T titanium 1x2 teeth 0.5 mm

- Straight shafts
- Serrated handle
- Polished finish
- Overall length 105 mm
- Available in Titanium or Stainless Steel

# 10L FOLDERS & INSERTERS

## Folding Forceps For Acrylic & Silicone Lens



### F-2106T titanium *With Lock*

- Angled shafts
- Round handle
- Overall length 110 mm

## Steinert Folding Forceps For Acrylic & Silicone Lens



### F-2107T titanium

- Angled shafts
- Flat handle
- Overall length 110 mm

## Insertion Forceps For Acrylic Lenses



### F-2132T titanium

- Flat handle
- Overall length 102 mm

## Buratto II Inserter For Acrylic Lens



### F-2135T titanium

- Flat handle
- Cross action
- Overall length 107 mm

# 10L FORCEPS

## T- END Lens Forceps



### F-2010T titanium

- Curved shaft
- Round handle
- Polished finish
- Overall length 120 mm

## Curved IOL Holding Forceps



### F-2013T titanium

- Curved shaft
- Flat handle
- Polished finish
- Overall length 85 mm

## Dementiev PRL Implantation Forceps, Angled



### F-2015T titanium

Designed to implant PRL

- Round knurled handle
- Overall length 120 mm

## IOL Cartridge Loading Forceps



### F-2141T titanium

Designed to implant PRL

- Round knurled handle
- Overall length 120 mm

## IRIS FORCEPS

### Iris Forceps With Serration



#### F-100T titanium

- Flat serrated handle
- Fine serrated tips
- Overall length 72 mm

### Iris Forceps With Teeth



#### F-101T titanium

- Straight shaft
- 1x2 teeth
- Flat serrated handle
- Polished finish
- Overall length 72 mm

### Iris Forceps With Teeth



#### F-102T titanium

- Curved shaft
- 1x2 teeth
- Flat serrated handle
- Polished finish
- Overall length 72 mm

## MOSQUITO FORCEPS

### Halsted Straight Hemostatic Mosquito Forceps



#### F-122S stainless steel

- Straight serrated jaws
- Ring handle
- Dull finish
- Overall length 120 mm

### Halstead Curved Hemostatic Mosquito Forceps



#### F-123S stainless steel

- Curved shafts
- Pointed tips
- Tying platform
- Smooth handle
- Polished finish
- Overall length 110 mm

# JEWELER FORCEPS

**Jeweler Forceps, Straight**



**F-111T titanium**

- Straight shafts
- Pointed tips
- Tying platform
- Smooth handle
- Polished finish
- Overall length 110 mm

**Jeweler Forceps, Straight**



**F-113T titanium**

- Straight shafts
- Pointed tips
- Tying platform
- Smooth handle
- Polished finish
- Overall length 110 mm

**Jeweler Forceps, Curved**



**F-115T titanium**

- Curved shafts
- Pointed tips
- Tying platform
- Smooth handle
- Polished finish
- Overall length 110 mm

# SUPERIOR RECTUS FORCEPS

**Troutman Superior Rectus Forceps**



**F-137T titanium**

- Angled with 1x2, 0.5mm teeth
- Shafts angled, 10 mm from tip to bend
- Serrated handle
- Polished finish
- Overall length 107 mm

**Elsching Superior Rectus Forceps**



**F-138T titanium**

- Straight with 1x2, 0.5mm teeth
- Straight shafts
- Serrated handle
- Polished finish
- Overall length 107 mm

# SUTURING FORCEPS

## Castroviejo Suturing Forceps



- F-0600T 1x2 teeth 0.12mm
- F-0601T 1x2 teeth 0.3 mm
- F-0602T 1x2 teeth 0.5 mm
- F-0603T 1x2 teeth 1.0 mm

- 6 mm tying platforms
- Dull finish
- Flat serrated handle
- Overall length 108 mm

## Castroviejo Angled Suturing Forceps



### F-0605T titanium

- 1x2 teeth 0.12 mm
- Angled shaft
- 6 mm tying platforms
- Dull finish
- Flat serrated handle
- Overall length 108 mm

## Paufique Suturing Forceps



### F-0606T titanium

- 1x2 teeth 0.5 mm
- 6 mm tying platforms
- Flat serrated handle
- Dull finish
- Overall length 87 mm

## Suture Removing Forceps



### F-0625T titanium

- 1x2 teeth 0.5 mm
- 6 mm platforms
- Flat serrated handle
- Dull finish
- Overall length 90 mm

# TYING FORCEPS

## Kelman-McPherson Tying Forceps



F-090T Tying Platforms 4 mm  
 F-091T Tying Platforms 8 mm  
 F-093T Tying Platforms 6 mm

- Titanium
- Flat serrated handle
- Overall length 85 mm

## McPherson Straight Tying Forceps



F-171T titanium

- Straight shafts
- Tying platforms 4 mm
- Flat serrated handle
- Polished finish
- Overall length 85 mm

## Kelman-McPherson Angled Tying Forceps, Long



F-174T Tying Platforms 8 mm  
 F-175T Tying Platforms 10 mm

- Flat serrated handle
- Overall length 110 mm

## McPherson Curved Tying Forceps, Long Size



F-177T titanium

- Curved shafts
- Tying platforms 6 mm
- Flat serrated handle
- Polished finish
- Overall length 110 mm

## Catalano Curved Tying Forceps



F-182T titanium

- 6 mm V-groove tying platform
- Round serrated handle
- Overall length 108 mm

## Tennant Straight Tying Forceps



F-185T titanium

- Extra delicate smooth jaws
- 6 mm tying platform
- For 9-0 to 11-0 sutures
- Round serrated handle
- Overall length 108 mm

# HOOKS

## Iris Hook, Angled



## H-020 Titanium Handle

- 1 mm hook shaped tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

## Lewicky Hook



## H-021 Titanium Handle

- 0.2 mm hook on curved 10 mm shaft
- Round serrated handle
- Dull finish
- Overall length 130 mm

## Kuglen Hook



## H-030 Titanium Handle

- "H" shape tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

## Sinsky Hook



## H-032 Titanium Handle

- 10 mm from tip to bend
- Round serrated handle
- Dull finish
- Overall length 130 mm

## Lester Lens Manipulator



## H-0331 Titanium Handle

- Angled shaft with 0.2 mm hourglass shape tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

## Bechert Nucleus Rotator



## H-034 Titanium Handle

- "Y" shape tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

# HOOKS

## Rentsch Boat Hook



### H-035 Titanium Handle

- Round serrated handle
- Dull finish
- Overall length 130 mm

## Fenzl Hook



### H-036 Titanium Handle

- Round serrated handle
- Dull finish
- Overall length 130 mm

## Delicate Muscle Hook



### H-0405 10 mm

- Small bulbous tip
- Flat serrated handle
- Dull finish
- Overall length 140 mm

## Gass Retinal Detachment Hook



### H-060

- 13 mm flat hook
- 1.5 mm hole
- Serrated handle
- Dull finish
- Overall length 140 mm

## DIAMOND KNIVES


### UNIVERSAL STRAIGHT KNIVES

**K-10**

Round knurled handle  
Overall length 110 mm



**DK050**  
45° Single Edge  
1.00 mm width



**DK052**  
45° Double Edge  
1.00 mm width



**DK053**  
Trifacet  
1.00 mm width



### ANGLED PHACO KNIVES

**K-20**


Angled Titanium Handle  
Overall length 130 mm



**Trapezoid  
Self-Diving Blades**

DK101 2.7/3.2 mm width  
DK104 2.3/2.8 mm width


Profile of the incision  
directed to the Anterior  
Chamber



**Facet-Free  
Crescent Blade**


DK091 1.5 mm width

Facets-free blades facilita-  
tes very smooth controlled  
cut.



**Clear Cornea Self-  
Diving Blades**

DK072 2.8 mm width  
DK074 3.2 mm width




# SAPPHIRE KNIVES

## UNIVERSAL STRAIGHT KNIVES


**K-10**  
 Round knurled handle  
 Overall length 110 mm




**SK050**  
 45° Single Edge  
 1.00 mm width



**SK052**  
 45° Double Edge  
 1.00 mm width



**SK053**  
 Trifacet  
 1.00 mm width



**SK070**  
 Lancet 60°  
 1.00 mm width



## ANGLED PHACO KNIVES


**K-20**  
 Angled Titanium Handle  
 Overall length 130 mm



### Trapezoid Blades


**SK301** 2.5/2.9 mm width  
 All Sharp

**SK302** 2.75/3.2 mm  
 width  
 All Sharp



### Crescent Blades


**SK092** 2.0 mm width



### Clear Cornea Blades

**SK701** 2.8 mm width

**SK703** 3.2 mm width



# PHACOEMULSIFICATION

## NUCLEUS CRACKERS / PRECHOPPERS

**Ernest Nucleus Cracker**



**P-025 Flat Handle**

- Cross action
- Overall length 102 mm

**Prechopper Angled**



**P-111 stainless steel**

- Cross action
- Round handle
- Maximum opening 2.8 mm
- Overall length 102 mm

**Prechopper Curved**



**P-112 stainless steel**

- Cross action
- Round handle
- Maximum opening 2.8 mm
- Overall length 102 mm

**Combo Prechopper Straight**



**P-116 stainless steel**

- Specially designed jaws to split soft and hard nucleus with sharp blades on the upper side of the tips and blunt on the other side
- Cross action
- Round handle
- Maximum opening 2.8 mm
- Overall length 102 mm

# PHACOEMULSIFICATION

## NUCLEUS DIVIDERS

### Kansas Nucleus Fragment Removing Forceps



#### P-0201 titanium

Designed to remove fragments of the nucleus through a small incision. Has two rows of delicate teeth on jaws

- Flat serrated handle
- Overall length 107 mm

### McIntyre Nucleus Bisector



#### P-060

- Round serrated handle
- Dull finish
- Overall length 130 mm

### McIntyre Nucleus Trisectort



#### P-061

- Round serrated handle
- Dull finish
- Overall length 130 mm

### McIntyre Nucleus Cutting Board



#### P-062

- Round serrated handle
- Dull finish
- Overall length 130 mm

# PHACOEMULSIFICATION

## PHACO CHOPPERS

### Nagahara Phaco Chopper



**P-063 Right Handed Surgeon**

**P-064 Left Handed Surgeon**

- Round serrated handle
- Dull finish
- Overall length 130 mm

### Rosen Phaco Chopper Universal



**P-065**

- Round serrated handle
- Dull finish
- Overall length 130 mm

### Phaco Cleaver



**P-068 Right Handed Surgeon**

**P-069 Left Handed Surgeon**

- Round serrated handle
- Dull finish
- Overall length 130 mm

### Quick Chopper



**P-125**

- Designed to quickly and easily chop all grades of nuclei
- 0.75 mm tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

# PHACOEMULSIFICATION

## NUCLEUS MANIPULATORS

### Koch Nucleus Spatula



#### P-070

- Gently curved
- Duck bill tip with notches
- Round serrated handle
- Overall length 130 mm

### Small Pupil Snapper Hook & Phaco Chopper



#### P-077

- Round serrated handle
- Dull finish
- Overall length 130 mm

### Rowen Phaco Spatula (Rosen Phaco Splitter & Bechert Nucleus Rotator)



#### P-095

- Round serrated handle
- Dull finish
- Overall length 145 mm

### Nucleus Sustainer



#### P-096

- Second hand instrument for hard cataracts, to be used with pre-choppers
- Round serrated handle
- Dull finish
- Overall length 130 mm

### Micro Finger



#### P-097

- Round serrated handle
- Dull finish
- Overall length 130 mm

# PHACOEMULSIFICATION

## MISCELLANEOUS

### Irrigation Handpiece For Bimanual Technique



**P-081**  
 21 GA  
 2x0.35 mm side ports  
 Regular, curved tip,  
 w/out polisher

### Aspiration Handpiece For Bimanual Technique



**P-082**  
 21 GA  
 1x0.35 mm top port  
 Regular, curved tip,  
 w/out polisher

### Aspiration Handpiece For Bimanual Technique



**P-0821**  
 22 GA  
 1x0.35 mm top port  
 Polisher, curved tip

### Phaco Chopper Rosen Type With Irrigation



**P-065/I**

- Irrigating handle
- 20 G tip with Rosen type phaco chopper
- Two 0.5mm irrigating ports
- Round handle
- Dull finish
- Overall length 120 mm

### Anterior/Posterior Capsule Polisher



**7-101**

- Angled 135 degree
- Round serrated handle
- Dull finish
- Overall length 130 mm

# NEEDLE HOLDERS

**Barraquer Needle Holder, Standard Jaws 12 mm**



**N-013 *Straight w/out Lock***

- Overall length 100 mm
- Round serrated handle

**Barraquer Needle Holder, Extra Fine Jaws 8mm**



**N-025 *Curved w/out Lock***

- Overall length 100 mm
- Round serrated handle

**Barraquer Needle Holder, Standard Jaws 12 mm**



**N-031 *Curved w/out Lock***

- Overall length 115 mm
- Round serrated handle


**Barraquer Needle Holder, Fine Jaws 12 mm**



**N-040 *Curved with Lock*  
N-041 *Curved w/out Lock***

- Overall length 115 mm
- Round serrated handle

**Barraquer Needle Holder, Extra Fine Jaws 8mm**



**N-045 *Curved with Lock***

- Overall length 115 mm
- Round serrated handle

**Barraquer Needle Holder, Strong Jaws 12 mm**



**N-090 *Curved with Lock*  
N-091 *Curved w/out Lock***

- Titanium
- Overall length 125 mm
- Round knurled handle

## NEEDLE HOLDERS

**Kalt Needle Holder Standard, Curved**



**N-081**

- Titanium
- 12 mm curved jaws
- Flat serrated locking handle
- Polished finish
- Overall length 140 mm

## PROBES & DILATORS

**Bowman Lachrymal Probe**



**D-010 size 0000 & 000**

**D-011 size 00 & 0**

**D-012 size 1 & 2**

**D-013 size 3 & 4**

**D-014 size 5 & 6**

**D-015 size 7 & 8**

**Lachrymal Pigtail Probe**



**D-030**

- 12 mm pigtail curved probes with eyes
- Round serrated handle
- Polished finish
- Overall length 140 mm

**Wilder Lachrymal Dilator**



**D-050 size 1, 19mm**

**D-051 size 2, 23mm**

**D-052 size 3, 32mm**

- Blunt tip
- Round knurled handle
- Polished finish
- Overall length 100 mm

**Castroviejo Double Ended Lachrymal Dilator**



**D-060 size 1 & 2**

- Round knurled handle
- Polished finish
- Overall length 140 mm

**Muldoon Lachrymal Dilator**



**D-061**

- Medium 27 mm taper distal end is straight to act as a guide for the tapered dilator
- Round knurled handle
- Polished finish
- Overall length 100 mm

# RETRACTORS

## Iris Retractor, Right



### R-060

- Round serrated handle
- Polished finish
- Overall length 100 mm

## Iris Retractor, Left



### R-061

- Round serrated handle
- Polished finish
- Overall length 100 mm

## Small Pupil Iris Retractor



### R-0701

*Plastic  
Disposable,  
Set of 4*

## Knapp Lachrymal Sac Retractor



### R-014

- Blunt 4 mm x 8 mm four prong tip
- Flat serrated handle
- Polished finish
- Overall length 130 mm

## Desmarres Lid Retractor



### R-020 size 0 1.0 mm

R-021 size 1 13.0 mm

R-022 size 2 15.0 mm

R-023 size 3 17.0 mm

- Round serrated handle
- Polished finish
- Overall length 130 mm

# SCISSORS

## CORNEAL SCISSORS

**Castroviejo Curved Corneal Scissors, Left**



**S-010**

- Left
- 11 mm blades
- Blunt tips
- Flat handle
- Polished finish
- Overall length 100 mm

**Castroviejo Curved Corneal Scissors, Right**



**S-0101**

- Right
- 11 mm blades
- Blunt tips
- Flat handle
- Polished finish
- Overall length 100 mm

**Castroviejo Universal Corneal Scissors**



**S-011 7.5 mm Blades**

**S-012 11 mm Blades**

- Flat Handle

**S-034 7.5 mm Blades**

**S-035 11 mm Blades**

- Round handle
- Blunt tips
- Polished finish
- Overall length 110 mm

**Katzin Curved Corneal Transplant Scissors, Left**



**S-020**

- Left
- Strongly curved
- Blunt tips
- Flat handle
- Polished finish
- Overall length 100 mm

**Katzin Curved Corneal Transplant Scissors, Right**



**S-0201**

- Right
- Strongly curved
- Blunt tips
- Flat handle
- Polished finish
- Overall length 100 mm

# SCISSORS

## CORNEAL SCISSORS

**Castroviejo Corneal Section Scissors, Right**



**S-024**

- Lower blade 0.5 mm longer than upper blade
- 14 mm blades
- Flat handle
- Polished finish
- Overall length 106 mm

**Castroviejo Corneal Section Scissors, Left**



**S-0241**

- Lower blade 0.5 mm longer than upper blade
- 14 mm blades
- Flat handle
- Polished finish
- Overall length 106 mm

## STITCH SCISSORS

**Westcott Stitch Scissors**



**S-046 Flat Handle**

- 13 mm blades
- Sharp pointed tips
- Polished finish
- Overall length 115 mm

**Westcott Stitch Scissors**



**S-047 Round Handle**

- 13 mm blades
- Sharp pointed tips
- Polished finish
- Overall length 115 mm

## TENOTOMY SCISSORS

**Westcott Curved Tenotomy Scissors**



**S-042 Flat Handle**

- 16 mm Blades
- Blunt tips
- Polished finish
- Overall length 120 mm

# SCISSORS

## CAPSULAR SCISSORS

**Side Port Capsulotomy Scissors, Left**



**S-0374 Left**

- Intraocular titanium handle
- 20 Gauge

**Side Port Capsulotomy Scissors, Right**



**S-0375 Right**

- Intraocular titanium handle
- 20 Gauge

**Clayman-Vannas Swan Neck Straight Scissors**



**S-0501**

- Sharp pointed tips
- 7 mm from mid-screw to tip
- Smooth handle
- Polished finish
- Overall length 83 mm

**Vannas Capsulotomy Scissors**



**S-050 Flat Handle**

- 6 mm blades
- Sharp pointed tips
- Polished finish
- Overall length 85 mm

**Vannas Curved Capsulotomy Scissors**



**S-052**

- Sharp pointed tips
- 7 mm from mid-screw to tip
- Smooth handle
- Polished finish
- Overall length 83 mm

**Gills-Vannas Straight Capsulotomy Scissors**



**S-056**

- 10 mm blades
- Sharp pointed tips
- Flat handle
- Polished finish
- Overall length 85 mm

# SCISSORS

## CAPSULAR SCISSORS

**Gills-Vannas Side Angled Capsulotomy Scissors**



**S-057**

- 10 mm blades
- Sharp pointed tips
- Flat handle
- Polished finish
- Overall length 85 mm

**Gills-Vannas Curved Capsulotomy Scissors**



**S-058**

- 10 mm blades
- Sharp pointed tips
- Flat handle
- Polished finish
- Overall length 85 mm

**Gills-Vannas Angled Capsulotomy Scissors**



**S-0581**

- 10 mm blades
- Sharp pointed tips
- Flat handle
- Polished finish
- Overall length 85 mm

**McPherson-Vannas Straight Iris Scissors**



**S-062**

- 8 mm blades
- Sharp pointed tips
- Round handle
- Polished finish
- Overall length 85 mm

## SCISSORS

### RING HANDLE SCISSORS

#### Straight Eye Scissors



#### S-080

- Pointed tips 28 mm
- Ring handle
- Polished finish
- Overall length 125 mm

#### Curved Eye Scissors



#### S-081

- Pointed tips 28 mm
- Ring handle
- Polished finish
- Overall length 125 mm

#### Curved Enucleation Scissors



#### S-091

- Curved round blades 48 mm
- Blunt tips
- Ring handle
- Polished finish
- Overall length 148 mm

# VITREORETINAL INSTRUMENTS

## VERTICAL SCISSORS



**V-202/V-003**  
Vertical 70°  
Long Sharp Blades

For pilling and cutting of fine epiretinal membranes in plane perpendicular to the retina

## STRAIGHT SCISSORS



**V-211/V-003**  
Straight  
Rounded Blades

Rounded blades allow safer cutting near the retinal surface

## HORIZONTAL SCISSORS 55°



**V-206/V-001**  
Horizontal 55°  
Medium 1.7 mm Blades

For cutting in the plane parallel to the retina and cutting membranes adherent to the retinal surface

## HORIZONTAL SCISSORS 25°



**V-212/V-001**  
Horizontal 25°  
Medium 1.7 mm Blades

For cutting in the plane parallel to the retina and cutting membranes adherent to the retinal surface

## ANGLED SCISSORS



**V-208/V-003**  
Angled 45° Horizontal  
Regular Blades

**V-2084/V-003**  
Angled 45° Horizontal  
Regular Blades With  
Illumination

## DIAMONIZED STRAIGHT GRIPPING FORCEPS



**V-301/V-003**  
Regular Size, 20GA

**V-3012/V-003**  
Regular Size, 23GA

**V-3013/V-003**  
Long Size, 20GA

# VITREORETINAL INSTRUMENTS

## POINTED MICROFORCEPS WITH PLATFORM



**V-327/V-003**  
Pointed End Gripping  
Forceps, 20GA

For manipulation of the retina

## DIAMONIZED ANGLED GRIPPING FORCEPS



**V-303/V-003**  
Regular Size, 20GA

**V-3032/V-003**  
Regular Size, 20/23GA

For membrane peeling and strip-  
ping

## HOCKEY STICK GRIPPING FORCEPS



**V-309/V-003**  
Regular Size, 20GA

For membrane peeling and strip-  
ping

## "PIC" FORCEPS



**V-3251/V-003**  
Long Size, 23GA

For grasping fine membranes  
adherent to the retinal surface

## FOREIGN BODY FORCEPS



**V-321/V-001**  
Spring Gripping, 20GA

For foreign body removing

## SUBRETINAL FORCEPS



**V-210/V-003**  
Angled

**V-2101/V-003**  
Angled  
With Illumination

# VITREORETINAL INSTRUMENTS

## SUBRETINAL FORCEPS, THOMAS TYPE



**V-403/V-003**  
Angled

For peeling and stripping of subretinal membranes. Useful for manipulation in small folds (macular puckers).


## HORIZONTAL SUBRETINAL FORCEPS



**V-407/V-003**  
Horizontal  
4.5mm Tips

Subretinal forceps are used for grasping membranes through retinotomy

## Silicone Oil Infusion Cannula



**V-025 4.0mm**  
**V-026 6.0mm**

- Stainless steel 20 G infusion cannulas with suture flange
- Each cannula is complete with silicone tubing and luer lock

## Backflush Instrument With Active Aspiration




**V-600**


- Overall length 130 mm

**V-600 T**

Silicone tubing for backflush instrument with active aspiration, 5pcs/box



**V-610**  
Curved needle



**V-611**  
Hooked needle



**V-618**  
Silicone tipped  
needle curved



**V-619**  
Silicone brush  
needle curved

# VITREORETINAL INSTRUMENTS

## Flushing Adapter For Vitreo Instruments



### V-000

- Titanium
- To flush vitreoretinal instruments with distilled water after each surgery

## Vitreoretinal Instrument Handle, Std. Model



### V-001T

Universal model with adjustable mechanism inside

- Round knurled titanium handle
- One finger control
- Overall length 110 mm

## Vitreoretinal Instrument Handle, Squeeze Model



### V-003T Angled

- Round titanium handle
- Two finger control
- Overall length 120 mm

# VITREORETINAL SPATULAS

## Membrane Spatula 2.5mm



### I-084

- Round handle
- Overall length 130 mm

## Membrane Spatula Knife



### I-090

- Round handle
- Overall length 130 mm

## Delicate Membrane Pick, 20 G




### I-099

- 0.9 mm membrane pick angled at 120° and tapers to 0.3mm
- Round handle
- Overall length 130 mm

# SPATULAS

## Straight Spatula



I-010 width 0.25mm, length 10mm  
 I-011 width 0.4mm, length 14mm  
 I-013 width 0.8mm, length 14mm  
 I-015 width 1.5mm, length 12mm  
 I-016 width 2.0mm, length 10mm

- Overall length 120 mm


## Curved Spatula



I-020  
*width 0.25 mm, length 12mm*  
 I-021  
*width 0.40 mm, length 12mm*

- Round handle
- Overall length 120 mm

## Double Spatula



I-031 *width 0.5mm and 1.0mm, length 12mm*  
 I-032 *width 0.30mm and 0.5mm, length 12mm*

- Round handle
- Overall length 130 mm

## Castroviejo Double Ended Cyclodialysis Spatula



I-050  
*0.5 mm by 10 mm and 0.5 mm by 15 mm spatulas*

- Round handle
- Overall length 130 mm

## Castroviejo Cyclodialysis Spatula



I-051  
*1 mm by 10 mm spatula*  
 I-052  
*0.75 mm by 10 mm spatula*

- Round handle
- Overall length 130 mm

## Dementiev PRL Microspatula



## I-140

Designed to place the haptics of the PRL under the iris  
 • 12x0.60 mm, very fine  
 • Tip covered with tungsten carbide

- Round serrated handle
- Overall length 120 mm

# SPATULAS

## LASEK/PRK Hockey Knife



### I-001

- Round handle
- Overall length 130 mm

Hockey knife can be used to remove epithelium for PRK procedure and during LASIK retreatment

## LASIK Flap Spatula



### I-005

- Round handle
- Overall length 130 mm

## Paton Spatula & Spoon




### I-110

- Round handle
- Overall length 145 mm

Spatula can be used to remove epithelium for PRK procedure

# SPECULUM


## Temporal Speculum, Lieberman Style



**L-040 Adult Size**

*Round Branches  
Titanium  
"V" Shape Blades*

## Temporal LASIK Speculum, Lieberman Style



**L-040 L**

*Round Branches  
Titanium  
Specially designed for pivoting  
microkeratome with nasal position  
of the microkeratome head  
Rounded Blades*

## Nasal Speculum, Lieberman Style



**L-0411 Adult Size**

*Flat Branches  
Titanium  
"V" Shape Blades*

## Nasal Speculum, Lieberman Style



**L-0411K**

*Flat Branches  
Titanium  
Straight Blades*


## Nasal LASIK Speculum, Lieberman Style



**L-0411L**

*Flat Branches  
Titanium  
Specially designed for pivoting  
microkeratome with temporal posi-  
tion of the microkeratome head  
Rounded Blades*

## Barraquer Eye Speculum




**L-022 Adult, 14mm Blades  
L-023 Child, 11mm Blades**

- Blade spread 18 mm
- Overall length 40 mm

# SPECULUM


## Barraquer Eye Speculum Solid Blades



### L-0221 Adult 14 mm Blades

- Blade spread 18 mm
- Overall length 40 mm

## Barraquer Eye Speculum Open Blades



### L-025 Adult Size

### L-026 Child Size

- Blade spread 18 mm
- Dull finish
- Overall length 40 mm

## Sauer Lid Speculum, Child



### L-030

- Solid 8 mm blades
- Blade spread 20 mm
- Overall length 35 mm

## Sauer Lid Speculum, New Born



### L-031

- Solid 5 mm blades
- Blade spread 20 mm
- Overall length 25 mm


## Sauer Lid Speculum, Premature



### L-032

- Solid 4 mm blades
- Blade spread 20 mm
- Overall length 20 mm

## Castroviejo Speculum



### L-020

- 15 mm fenestrated blades
- Screw controlled spread of blades
- Blade spread 30 mm
- Dull finish
- Overall length 80 mm

# SPECULUM


**Temporal Speculum With Aspiration, Lieberman Style**



**L-080 A Adult Size**  
**14mm "V" Shape Blades**

Six ports on each blade allow continuous aspiration for a debris free stroma resulting in better postoperative vision


**Temporal Speculum With Aspiration, Lieberman Style**



**L-082 A Child Size**  
**11mm "V" Shape Blades**

Six ports on each blade allow continuous aspiration for a debris free stroma resulting in better postoperative vision

**Temporal Speculum With Aspiration, Lieberman Style**



**L-080L A**  
**Specially designed for LASIK**  
**Rounded Blades**

Six ports on each blade allow continuous aspiration for a debris free stroma resulting in better postoperative vision

# CANNULAS

## Anterior Chamber (AC) Cannula (Rycroft)



**A-051-23 23G**  
Angled 45 degree  
Overall length  
22 mm

## AC Washout Cannula



**A-055-16 16G**  
30 degree angled,  
Spatulated, 6mm  
from end length  
25mm

## Anterior Chamber Cannula (Knolle)




**A-063-23 23G**  
2 mm angled tip,  
length 23mm

## Anterior Chamber Maintainer (Lewicky)



**A-065**  
Silicone tubing  
attached to 20 G,  
3.5mm cannula,  
self-retained in AC

## Anterior Chamber Maintainer (Girard)




**A-067**  
Silicone tubing at-  
tached to 23G, 5mm  
cannula with Atkin-  
son point, slightly  
roughened for self-  
retaining in A C

## Buratto LASIK Irrigating Cannula



**A-379-23 23G**  
One front and two  
side ports.  
Designed for equal  
flow in all direc-  
tions.

## Cortex Extractor Cannula (Simcoe)




**A-091-21-21G**  
**A-091-13-21G**  
Smooth closed tip  
with top aspirating  
port into curve open-  
ing

## Cortex Extractor Cannula



**A-1098**  
For Irrigation or  
aspiration, 25G

## I & A Cannula (Simcoe)



**A-129-0.3 23/23G 0.3 mm port**  
**A-129-0.4 23/22G 0.4 mm port**

# CANNULAS

**Silicone Bulb**

**Adapter**



**A-301/A-303**

**Silicone bulb with adapter**

**Microsurgical Handle**



**A-307**

**Standard Male / Male  
Luer connector**

**Microsurgical Handle**



**A-308**

**Standard Male / Female Luer  
connector**

# MISCELLANEOUS

- M-0300 6 mm
- M-0301 6,5 mm
- M-0302 6,75 mm
- M-0303 7 mm
- M-0304 7,25 mm
- M-0305 7,5 mm
- M-0306 7,75 mm
- M-0307 8 mm
- M-0308 8,25 mm
- M-0309 8,5 mm
- M-0310 9 mm
- M-0311 9,5 mm

## Corneal Trephine



## M-1320

Consists of Universal Adjustable Trephine Handle and Trephine Blades of various sizes

Please, specify blade size when ordering

## Universal Handle



## Thermocautery



## M-000

- 6 mm ball tip
- Titanium handle

## Kelly Decement Membrane Punch



## M-011

- Produce clean 0.75 mm round hole without tissue tags
- Serrated squeeze action handle
- Titanium handle
- Stainless steel tip

## Enucleation Spoon



## M-060

- Round serrated handle
- Overall length 145 mm

## Evisceration Spoon



## M-062 large size

- Round knurled handle
- Overall length 130 mm

## Meyerhoeffer Chalazion Curette



- M-063 size 0 - 1.5 mm
- M-064 size 1 - 1.75 mm
- M-065 size 2 - 2.0 mm
- M-066 size 3 - 2.5 mm
- M-067 size 4 - 3.50 mm

# MISCELLANEOUS

**Veber Modified Lens Loop**



**M-071**

- 6 mm x 9 mm loop with cross hair
- Round knurled handle
- Overall length 130 mm

**Wilder Lens Loop**



**M-072**

- Round knurled handle
- Overall length 130 mm

**Lid Plate**



**M-50**

- 20 mm and 24 mm wide
- Polished finish
- Overall length 110 mm

**Fine - Thornton Phaco Fixation Ring**



**M-0341 14 mm diameter ring with teeth**

**M-036 13 mm diameter ring with teeth**

- 3/4 open ring with swivel
- Round handle
- Overall length 90 mm


**Towel Clamp,**



**M-081**

- Ring handle
- Overall length 90 mm

**Two Handed Screw Injector for Acrylic IOLs**



**M-281**

Designed to insert foldable IOL into the cartridge

## BASIC EYE SET

#	ITEM	DESCRIPTION	Q-TY
1	B-010	Blade Holder	1
2	C-010	Castroviejo Caliper	1
3	F-100	Straight Iris Forceps with Serration	1
	F-101	Straight Iris Forceps with Teeth	
	F-102	Curved Iris Forceps with Teeth	
4	F-050	Colibri Corneal Forceps	1
	F-0501	Colibri Corneal Forceps	
	F-053	Colibri Corneal Forceps	
	F-054	Colibri Corneal Forceps	
	F-0541	Colibri Corneal Forceps	
5	F-0822	Castroviejo Fixation Forceps	1
	F-0831	Castroviejo Fixation Forceps	
6	F-058	Bonn Corneal Forceps	1
	F-055	Catalano Corneal Forceps	
	F-0600	Castroviejo Suturing Forceps 0.12mm	
	F-0605	Castroviejo Suturing Forceps 0.12mm	
7	F-0601	Castroviejo Suturing Forceps 0.3mm	1
	F-0607	Bishop-Harmon Suturing Forceps	
8	F-0602	Castroviejo Suturing Forceps 0.5mm	1
	F-0606	Paufique Suturing Forceps 0.5mm	
9	F-070	Dressing Forceps	1
	F-071	Dressing Forceps	
	F-073	Bishop-Harmon Dressing Forceps	
10	F-2301	Fechtner Conjunctiva Forceps	1
	F-2302	Fechtner Conjunctiva Forceps	
11	F-171	McPherson Tying Forceps Straight	1
	F-177	McPherson Tying Forceps curved	
	F-185	Tennant Tying Forceps Straight	
	F-186	Tennant Tying Forceps Curved	
12	F-091	Kelman-McPherson Tying Forceps Angled	1
	F-092	Kelman-McPherson Tying Forceps Angled	
	F-174	Kelman-McPherson Tying Forceps Angled	
	F-175	Kelman-McPherson Tying Forceps Angled	
13	K-10/DK050	45° Single Edge Universal Diamond Knife	1
14	N-031	Barraquer Needle Holder Standard Curved	1
15	N-045	Barraquer Needle Holder Extra Fine Curved	1
16	R-020	Desmarres Lid Retractor size 0	1
17	R-021	Desmarres Lid Retractor size 1	1
18	R-022	Desmarres Lid Retractor size 2	1
19	R-23	Desmarres Lid Retractor size 3	1
20	S-042	Westcott Tenotomy Scissors Curved	1
	S-043	Westcott Tenotomy Scissors Curved	
21	S-045	Westcott Stitch Scissors	1
	S-047	Westcott Stitch Scissors	
22	S-0581	Gills-Vannas Capsulotomy Scissors Angled	1
23	L-040	Temporal Speculum Lieberman Style	1
24	L-022	Barraquer Eye Speculum	1
25	L-023	Barraquer Eye Speculum, Child Size	1
26	A-051	Cannula Rycroft	1
27	A-301	Adapter	1
28	A-303	Silicone Bulb	1
29	M-400	Titanium Sterilizing Tray with Silicone Finger Tip Mat	2

## CATARACT SET (ECCE)

#	ITEM	DESCRIPTION	Q-TY
1	B-010	Blade Holder	1
2	C-010	Castroviejo Caliper	1
3	F-050	Colibri Corneal Forceps	1
	F-0501	Colibri Corneal Forceps	
	F-053	Colibri Corneal Forceps	
	F-054	Colibri Corneal Forceps	
4	F-073	Bishop-Harmon Dressing Forceps	1
	F-070	Straight Dressing Forceps	
	F-071	Straight Dressing Forceps	
	F-072	Curved Dressing Forceps	
5	F-136	Superior Rectus Forceps	1
	F-137	Troutman Superior Rectus Forceps	
6	F-2001	Clayman IOL Forceps	1
	F-2002	Blades Angled Lens Forceps	
	F-2005	Lens Forceps	
	F-2006	Lens Holding Forceps	
	F-2007	Kratz Lens Holding Forceps	
	F-2009	Ring Lens Forceps	
	F-2010	T-End Forceps	
7	F-0600	Castroviejo Suturing Forceps 0.12mm	1
	F-0605	Castroviejo Angled Suturing Forceps 0.12mm	
	F-058	Bonn Corneal Forceps 0.12mm	
	F-055	Catalano Corneal Forceps 0.12mm	
8	F-0601	Castroviejo Suturing Forceps 0.3mm	1
	F-0607	Bishop-Harmon Suturing Forceps	
9	F-171	McPherson Tying Forceps Straight	1
	F-177	McPherson Tying Forceps Curved	
	F-185	Tennant Tying Forceps Straight	
	F-186	Tennant Tying Forceps Curved	
	F-1851	Tennant Tying Forceps Angled	
10	F-091	Kelman-McPherson Tying Forceps Angled	1
	F-092	Kelman-McPherson Tying Forceps Angled	
	F-174	Kelman-McPherson Tying Forceps Angled, Long	
	F-175	Kelman-McPherson Tying Forceps Angled, Long	
11	F-0301	Utrata Capsulorrhexis Forceps Straight	1
	F-0311	Utrata Capsulorrhexis Forceps Straight	
	F-0321	Utrata Capsulorrhexis Forceps Curved	
	F-0331	Utrata Capsulorrhexis Forceps Curved	
12	H-032	Sinsky Hook	1
13	H-030	Kuglen Hook	1
14	K-10/DK050	45° Single Edge Universal Diamond Knife	1
15	N-031	Barraquer Needle Holder Standard Curved	1
16	N-045	Barraquer Needle Holder Extra Fine Curved	1
17	S-010	Castroviejo Corneal Scissors Left	1
18	S-0101	Castroviejo Corneal Scissors Right	1
19	S-011	Castroviejo Universal Corneal Scissors	1
	S-012	Castroviejo Universal Corneal Scissors Long Blades	
	S-034	Castroviejo Universal Corneal Scissors	
20	S-042	Westcott Tenotomy Scissors Curved	1
	S-043	Westcott Tenotomy Scissors Curved	
21	S-0581	Gills-Vannas Angled Capsulotomy Scissors	1
22	S-045	Westcott Stitch Scissors	1
	S-047	Westcott Stitch Scissors	
23	I-031	Double Ended Spatula 0.5mm/1.0mm	1
	I-014	Straight Spatula 1.0mm	
24	L-022	Barraquer Wire Speculum	1
	L-0221	Barraquer Wire Speculum	
	L-025	Barraquer Wire Speculum	
25	L-040	Temporal Speculum Lieberman Style	1
	L-0411	Nasal Speculum Lieberman Style	
	L-0411K	Nasal Lieberman Style	
26	L-080A	Temporal Speculum Lieberman Style with Aspiration	1
27	15-129	I/A Cannula Simcoe	1
28	15-091	Cortex Extractor Cannula	1
	15-1098	Cortex Extractor Cannula	
29	15-301	Adapter	1
30	15-303	Silicone Bulb	1
31	M-072	Wilder Lens Loop	1
32	M-400	Titanium Sterilizing Tray with Silicone Finger Tip Mat	2

## PHACOEMULSIFICATION SET

#	ITEM	DESCRIPTION	Q-TY
1	C-010	Castroviejo Caliper	1
2	F-050	Colibri Corneal Forceps	1
3	F-073	Bishop-Harmon Dressing Forceps	1
4	F-136	Osher Superior Rectus Forceps	1
5	F-2001	Clayman IOL Forceps	1
6	F-0600	Castroviejo Suturing Forceps 0.12mm	1
7	F-0601	Castroviejo Suturing Forceps 0.3mm	1
8	F-182	Tennant Tying Forceps Curved	1
9	F-185	Tennant Tying Forceps Straight	1
10	F-0301	Utrata Capsulorhexis Forceps Straight	1
11	H-032	Sinsky Hook	1
12	H-034	Bechert Nucleus Rotator	1
13	H-030	Kuglen Hook	1
14	K-10/DK053	Trifacet Universal Straight Knife	1
15	K-20/DK092	Round Crescent Facet-Free Knife, Angled	1
16	K-20/DK071	2.5mm Clear Cornea Knife, Angled	1
	K-20/DK072	2.8mm Clear Cornea Knife, Angled	
	K-20/DK074	3.2mm Clear Cornea Knife, Angled	
17	K-20/DK101	2.7/3.2mm Trapezoid Self-Diving Diamond Knife, Angled	1
	K-20/DK103	2.5/2.9mm Trapezoid Self-Diving Diamond Knife, Angled	
	K-20/DK104	2.3/2.8mm Trapezoid Self-Diving Diamond Knife, Angled	
18	P-063	Nagahara Phaco-chopper, RHD	1
	P-064	Nagahara Phaco-chopper, LHD	
19	P-025	Ernest Nucleus Cracker	1
20	N-031	Barraquer Needle Holder Standard Curved	1
21	N-045	Barraquer Needle Holder Extra Fine Curved	1
22	S-042	Westcott Tenotomy Scissors Curved	1
23	S-045	Westcott Stitch Scissors	1
24	S-0581	Gills-Vannas Capsulotomy Scissors Angled	1
25	I-050	Castroviejo Cyclodialysis Double Ended Spatula	1
26	L-022	Barraque Wire Speculum	1
27	L-080A	Temporal Speculum Lieberman Style with Aspiration	1
28	A-129	I/A Cannula Simcoe	1
29	A-091	Cortex Extractor Cannula	1
30	A-301	Adapter	1
31	A-303	Silicone Bulb	1
32	P-081	Irrigation Handpiece for Bimanual Technique	1
33	P-082	Aspiration Handpiece for Bimanual Technique	1
34	P-110	Pre-Chopper, Straight	1
	P-111	Pre-Chopper, Angled	
	P-112	Pre-Chopper, Curved	
35	M-400	Titanium Sterilizing Tray with Silicone Finger Tip Mat	2

## CORNEAL TRANSPLANT SET

#	ITEM	DESCRIPTION	QTY
1	C-140	Corneal Transplant Marker	1
2	F-0814	Polack Corneal Transplant Fixation Forceps	1
3	F-0541	Castroviejo Colibri Corneal Forceps	1
4	F-0607	Bishop-Harmon Suturing Forceps	1
5	F-0600	Castroviejo Suturing Forceps 0.12mm	1
6	F-0601	Castroviejo Suturing Forceps 0.3mm	1
7	F-090	Kelman-McPherson Angled Tying Forceps	1
8	N-031	Barraquer Needle Holder Standard	1
9	N-045	Barraquer Needle Holder Extra Fine Curved	1
10	S-020	Katzin Curved Corneal Transplant Scissors Left	1
11	S-0201	Katzin Curved Corneal Transplant Scissors Right	1
12	S-024	Castroviejo Corneal Section Scissors Left	1
13	S-0241	Castroviejo Corneal Section Scissors Right	1
14	S-042	Westcott Tenotomy Scissors Curved	1
15	S-0581	Gills-Vannas Capsulotomy Scissors Angled	1
16	S-045	Westcott Stitch Scissors	1
17	L-022	Barraquer Wire Speculum	1
18	A-051	Cannula Rycroft	1
19	A-301	Adapter	1
20	A-303	Silicone Bulb	1
21	M-0305	Corneal Trephine Blade 7.5mm	1
22	M-0307	Corneal Trephine Blade 8.0mm	1
23	M-1320	Corneal Trephine Universal Handle	1
24	M-400	Titanium Sterilizing Tray With Silicone Finger Tip Mat	2

## LID SURGERY SET

#	ITEM	DESCRIPTION	Q-TY
1	B-010	Blade Holder	1
2	C-010	Castroviejo Caliper	1
3	F-0822	Castroviejo Fixation Forceps	1
4	F-1906	Desmarres Chalazion Forceps Oval, Large Size	1
5	F-1908	Lambert Chalazion Forceps Small Size	1
6	F-1912	Desmarres Chalazion Forceps	1
7	F-1901	Erhardt Lid Forceps	1
8	F-1920	Laser Lid Clamp Forceps, Large	1
9	F-1921	Laser Lid Clamp Forceps, Small	1
10	F-0600	Castroviejo Suturing Forceps, 0.12mm	1
11	F-0601	Castroviejo Suturing Forceps, 0.3mm	1
12	F-0602	Castroviejo Suturing Forceps, 0.5mm	1
13	F-073	Bishop-Harmon Dressing Forceps	1
14	F-2301	Fechtner Conjunctiva Forceps	1
15	F-136	Osher Superior Rectus Forceps	1
16	F-090	Kelman-McPherson Tying Forceps Angled	1
17	F-040	Ziegler Cilia Forceps	1
18	F-0740	Adson Dressing Forceps	1
19	F-0741	Adson Fixation Forceps	1
20	N-031	Barraquer Needle Holder Standard Curved	1
21	R-020	Desmarres Lid Retractor Size 0	1
22	R-021	Desmarres Lid Retractor Size 1	1
23	R-022	Desmarres Lid Retractor Size 2	1
24	R-023	Desmarres Lid Retractor Size 3	1
25	S-042	Westcott Tenotomy Scissors Curved	1
26	S-045	Westcott Stitch Scissors	1
27	A-051	Cannula Rycroft	1
28	A-301	Adapter	1
29	A-303	Silicone Bulb	1
30	M-50	Lid Plate	1
31	M-063	Meyerhoeffer Chalazion Curette Size 0, 1.5mm	1
32	M-064	Meyerhoeffer Chalazion Curette Size 1, 1.75mm	1
33	M-065	Meyerhoeffer Chalazion Curette Size 2, 2.0mm	1
34	M-066	Meyerhoeffer Chalazion Curette Size 3, 2.5mm	1
35	M-067	Meyerhoeffer Chalazion Curette Size 4, 3.5mm	1
36	M-400	Titanium Sterilizing Tray with Silicone Finger Tip Mat	2

## CHALAZION SET

#	ITEM	DESCRIPTION	Q-TY
1	B-010	Blade Holder	1
2	M-063	Meyerhoeffer Chalazion Curette Size 0, 1.5mm	1
3	M-064	Meyerhoeffer Chalazion Curette Size 1, 1.75mm	1
4	M-065	Meyerhoeffer Chalazion Curette Size 2, 2.0mm	1
5	M-066	Meyerhoeffer Chalazion Curette Size 3, 2.5mm	1
6	M-067	Meyerhoeffer Chalazion Curette Size 4, 3.5mm	1
7	F-1906	Desmarres Chalazion Forceps Oval Big Size	1
8	F-1908	Lambert Chalazion Forceps Small	1
9	F-1909	Lambert Chalazion Forceps Medium	1
10	F-1912	Lambert Chalazion Forceps	1
11	S-045	Westcott Stitch Scissors	1
12	M-400	Titanium Sterilizing Tray with Silicone Finger Tip Mat	1

### VITRECTOMY SET

#	ITEM	DESCRIPTION	Q-TY
1	C-010	Castroviejo Caliper	1
2	F-050	Colibri Corneal Forceps	1
3	F-073	Bishop-Harmon Dressing Forceps	1
4	F-0600	Castroviejo Suturing Forceps, 0.12mm	1
5	F-0601	Castroviejo Suturing Forceps, 0.3mm	1
6	F-171	McPherson Tying Forceps Straight	1
7	F-090	Kelman-McPherson Tying Forceps Angled	1
8	F-2202	Lambert Sleeve Spreading Forceps	1
9	F-2203	Landers Vitrectomy Lens Forceps	1
10	N-031	Barraquer Needle Holder Standard Curved	1
11	N-045	Barraquer Needle Holder Extra Fine Curved	1
12	S-0581	Gills-Vannas Angled Capsulotomy Scissors	1
13	S-042	Westcott Tenotomy Scissors Curved	1
14	S-045	Westcott Stitch Scissors	1
15	V-001	Vitreoretinal Instrument Handle, Standard Model	2
16	V-003	Vitreoretinal Instrument Handle, Squeeze Model	2
17	V-000	Flushing Adapter for Vitreoretinal Instruments	2
18	V-202	Vertical Scissors	1
19	V-211	Straight Scissors	1
20	V-206	Horizontal 55 Degree Scissors	1
21	V-301	Diamonized End-Gripping Forceps	1
22	V-303	Diamonized Angled Vitreoretinal Forceps	1
23	V-321	Foreign Body Forceps	1
24	V-600	Backflush Instruments With Active Aspiration	1
25	V-610	Curved Needle	1
26	V-611	Hooked Needle	
27	V-618	Silicone Tipped Needle, Curved	
28	V-619	Silicone Brush Needle, Curved	
24	I-061	Olk Angled Retinal Membrane Pick 110'	1
25	I-085	Membrane Spatula 5.0mm	1
26	I-090	Membrane Spatula Knife	1
27	L-040	Temporal Speculum Lieberman Style	1
28	M-400	Titanium Sterilizing Tray with Silicone Finger Tip Mat	2

### SUBRETINAL SET (ADD TO VITRECTOMY SET)

#	ITEM	DESCRIPTION	Q-TY
1	V-210	Subretinal Forceps	1
2	V-2101	Subretinal Forceps with Illumination	1

## NOTES

# CONTENTS

<b>Blade Holders, Calipers &amp; Markers</b>	<b>1</b>
<b>Forceps</b>	<b>2</b>
<b>Hooks</b>	<b>16</b>
<b>Knives</b>	<b>18</b>
<b>Phacoemulsification</b>	<b>20</b>
<b>Needle Holders</b>	<b>25</b>
<b>Probes &amp; Dilators</b>	<b>26</b>
<b>Retractors</b>	<b>27</b>
<b>Scissors</b>	<b>28</b>
<b>Vitreoretinal</b>	<b>33</b>
<b>Spatulas</b>	<b>37</b>
<b>Speculums</b>	<b>39</b>
<b>Cannulas</b>	<b>42</b>
<b>Miscellaneous</b>	<b>44</b>
<b>Instrument Sets</b>	<b>46</b>

[WWW.CILITA.COM](http://WWW.CILITA.COM)

**CILITA**

**45, VOKZALNAYA STR.,  
RYAZAN, RUSSIA, 390013  
Phone/Fax: 7 4912 983686  
E-MAIL: SALES@CILITA.COM**